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Instructions for use Secretion aspirator Allegra M30

Copyright © 2022 ASSKEA GmbH, Gebesee.

The safety of the **secretion aspirator Allegra M30** complies with the recognized rules of technology and meets the requirements of the **German Medical Devices Act**.

The **secretion aspirator Allegra M30** bears the **CE marking CE0494** in accordance with the EU Council Directive 93/42/EEC concerning medical devices and meets the essential requirements of Annex I of this directive and Article 120 of Regulation (EU) 2017/745 of the European Parliament and of the Council.

The **secretion aspirator Allegra M30** has been tested in accordance with IEC 62353.

The **quality management system** applied by ASSKEA GmbH is certified in compliance with the relevant international standards.

The **secretion aspirator Allegra M30** is a class IIa medical aspirator in accordance with EU Council Directive 93/42/EEC, Annex IX.

Errors and omissions excepted.

Contents

1	Information on delivery damage and customer feedback	47
2	User information	48
2.1	Using these instructions for use	48
2.2	Icons	48
2.2.1	Device, packaging and accessories.....	48
2.2.2	Operating elements	49
2.3	Symbol convention.....	49
2.4	Glossary	50
2.5	Purpose.....	51
2.5.1	Intended use	51
2.5.2	Essential functions.....	51
2.5.3	Applied parts	51
2.5.4	Indications	51
2.5.5	Contraindications	51
2.5.6	Precautionary measures	51
2.5.7	Restrictions on use	51
2.6	Basic safety instructions – ATTENTION!.....	52
2.7	User requirements	54
2.8	Information on product liability	54
3	Product description	55
3.1	Complete view of the ME system.....	55
3.1.1	secretion aspirator Allegra M30 with reusable secretion canister system.....	55
3.1.2	secretion aspirator Allegra M30 with disposable secretion canister system	55
3.2	Product contents	56
3.3	Product properties	56
3.4	Information on the bacterial filter and the adapter for disposable secretion canister system.....	57
3.5	Reusable secretion canister system for secretion aspirator Allegra M30	57
3.6	Disposable secretion canister system for secretion aspirator Allegra M30.....	58
3.7	Information on the battery.....	58
3.8	Warranty.....	59
4	Operation	60
4.1	Setup and startup.....	60
4.1.1	Connecting the secretion aspirator Allegra M30.....	60
4.2	Startup	61
4.2.1	Connecting the reusable secretion canister.....	61
4.2.2	Change of the bacterial filter or the adapter for the disposable secretion canister system	63
4.2.3	Connecting the disposable secretion canister system	64
4.3	Operation	67

4.3.1	Control elements and indicators of the secretion aspirator Allegra M30.....	67
4.3.2	Implementation of the aspiration.....	68
4.3.3	Replacement of the disposable liner „OneWay“	68
4.4	Decommissioning.....	69
4.4.1	secretion aspirator Allegra M30 with reusable secretion canister	69
4.4.2	secretion aspirator Allegra M30 with disposable secretion canister system	71
5	Maintenance	72
5.1	General information	72
5.2	Processing	73
5.2.1	Surface of the device.....	73
5.2.2	Reusable secretion canister system	73
5.2.3	External canister „Bag“	74
5.2.4	Holder for the external canister “Bag”	75
5.2.5	Suction tube	75
5.2.6	Adapter for the disposable secretion canister system	75
5.3	Further use of the device on a new patient	76
6	Troubleshooting	77
7	Transport, storage, and disposal	79
7.1	Transport / Return.....	79
7.2	Storage.....	79
7.3	Disposal	79
8	Technical data	80
9	Electromagnetic compatibility.....	81
10	Order information	82
11	Publishing information	83

1 Information on delivery damage and customer feedback

Caution!

Please check consignments for damage immediately upon delivery and remove all packaging. Any damage must be reported within three days; otherwise, these cannot be recognised.

The products described in this instruction for use are subject to constant development and improvement. For this reason, we welcome all customer feedback, comments and suggestions regarding our products and their accompanying documents, which contribute to improving the product, service or documentation.

en

Feel free to let us know what you think!

Use one of the following ways to let us know:

- ➔ by e-mail to info@loewensteinmedical.com
- ➔ by telephone at +49 (0) 2603/9600-0

2 User information

2.1 Using these instructions for use

Please read these instructions for use in their entirety before operating the **secretion aspirator Allegra M30** for the first time. The device is intended for operation by appropriately trained medical professionals, relatives of patients, as well as patients themselves provided that they, too, have been instructed accordingly. That being the case, all functions can be used safely. If support is required for the startup, operation or maintenance, please contact Löwenstein Medical (see section 11). Please also report unexpected operating performance or incidents to Löwenstein Medical (info@loewensteinmedical.com).



Please read the safety instructions (see section 2.6) to avoid hazards.

This instruction for use is a component of the **secretion aspirator Allegra M30**. Therefore please keep these instructions for use in an easily accessible place.

Include this instruction for use when passing the **secretion aspirator Allegra M30** on to third parties.

2.2 Icons

2.2.1 Device, packaging and accessories

Symbol	Meaning	Symbol	Meaning
	Attention: possible bodily injury, health risks or possible property damage		Note with useful information and tips
	Follow the instructions for use.		This device must not be disposed of in domestic waste.
	Power supply unit		Order number
	Protection class II		Serial number
	Protection class: Type BF (Body Floating)		Lot number
	Temperature limitation		Date of manufacture
	Air pressure limitation		Manufacturer
	Humidity limitation		Operating time (power-on time and power-off time)
	Do not use if packaging is damaged!		Do not reuse.
	International Protection: IP22 (see section 2.4)		Protect from moisture

Symbol	Meaning	Symbol	Meaning
	For reuse on one patient		Non-sterile
	Medical device		Unique Device Identification
	Distributor		Swiss Authorised Representative
	Importer		Follow the instructions for use

en

2.2.2 Operating elements

Symbol	Meaning	Symbol	Meaning
	Suction setting Low (vacuum ≤ -0.2 bar, flow rate 20 l/min)		LED status indicator for battery operation
	Suction setting Medium (-0.2 bar < vacuum ≤ -0.5 bar, flow rate 25 l/min)		LED status indicator for mains operation
	Suction setting High (-0.5 bar < vacuum ≤ -0.9 bar, flow rate 30 l/min)		On / Off button

2.3 Symbol convention

Symbol	Meaning
•	Bulleted list
1. 2.	Perform the actions in the specified order.

2.4 Glossary

A

Aspirate Aspirate is the generic term for secretions, bodily fluids and flushing liquids that are typically accumulated during aspiration of the upper airways. These can be aspirated easily using the device described here.

C

Contamination Contamination means that bacteria and viruses from the aspirate have come into contact with the device.

D

Degree of protection The degree of protection indicates the protection of applied parts against electric shock. Type BF applied parts have to be installed with an isolated ground and are not suitable for direct application to the heart.

I

incl. Abbreviation for "including"

IP22 Abbreviation for "International Protection / Protection Class"
The Protection Class specifies the degree of protection of the device against contact and ingress of liquids.

The **secretion aspirator Allegra M30** is protected against finger access and falling water drops at an inclination of up to 15°.

M

ME system Abbreviation for "medical electrical system"

MRI Abbreviation for "magnetic resonance imaging"
This technique generates sectional images of the human body with the aid of a very strong magnetic field for the purpose of analyzing the organs.

O

Overflow Overflow means that the aspirate is suctioned into the interior of the device.

P

Processing The processing procedure is required for each new patient. The term processing means that all parts that have (potentially) come into contact with the aspirate must be cleaned, disinfected and replaced if necessary. Reconditioning may be performed only by Löwenstein Medical or by an authorized service partner.

2.5 Purpose

2.5.1 Intended use

The **secretion aspirator Allegra M30** is a network-independent mobile medical aspirator and is used for temporary and primarily spontaneous aspiration of aspirate from the trachea.

Typical areas of application are:

- outpatient and inpatient care (professional healthcare facility environment),
- homecare (home healthcare environment), especially for the aspiration of aspirate from tracheostomy patients.

2.5.2 Essential functions

- Generation of vacuum (high vacuum)
- Generation of flow (high flow)

2.5.3 Applied parts

- There are no applied parts; however, the housing and suction tube are considered as type BF applied parts.

2.5.4 Indications

- Tracheostomy patients
- Aspiration due to impaired respiratory function
- Aspiration of blood, secretions, and food residues from the oral cavity, the pharyngeal zone and bronchial system

2.5.5 Contraindications

The **secretion aspirator Allegra M30** is contraindicated for the following applications:

- Liposuction
- Gynecology applications
- Dental applications
- Thoracic drainage
- Continuous drainage
- Applications in wound areas

2.5.6 Precautionary measures

Precautionary measures should be taken in the case of fresh oesophageal or tracheal anastomoses or tracheo-oesophageal fistulas.

2.5.7 Restrictions on use

- in medical rooms where potential equalisation is required (e.g. heart surgery)
- in potentially explosive areas / in the MRI environment
- outdoors / during transport

2.6 Basic safety instructions – ATTENTION!

Risk of damage due to improper power supply

Improper operation causes overvoltage in the device that may be transmitted to the operator.

- Ensure prior to startup that the mains supply for connecting the **secretion aspirator Allegra M30** is suitable at 100 V to 240 V AC and a mains frequency of 50-60 Hz.
 - Ensure prior to startup in UL listed markets such as the USA and Canada that the mains supply is designed for a supply voltage of 120 V AC.
 - Use only the supplied power supply unit with the **secretion aspirator Allegra M30** (see section 8).
-

Risk of application under non-approved conditions

- The device is intended for mobile, temporary outdoor use (according to IP protection class 22), but not for permanent outdoor use and operation outdoors.
-

Health risks due to handling of infectious or pathogenic germs

Infectious and pathogenic germs in the aspirate pose health risks.

- Never operate the device without a bacterial filter.
 - While carrying out cleaning / disinfection, make sure that the bacteria filter does not get wet.
 - Always perform the aspiration using a sterile disposable catheter.
 - The suction tube must never come into direct contact with the aspiration site.
 - Follow the hygiene, cleaning and decontamination instructions (section 5 Maintenance).
-

Risk of strangulation

- People may strangle themselves on the tubing or the mains cable, especially if tubes or cables are unduly long.
 - Ensure that no unauthorized / uninvolved persons are near the device during aspiration.
 - Store the device incl. accessories in the shipping carton until the next use.
-

Risk of damage due to electromagnetic phenomena

Medical electrical devices are subject to special precautionary measures regarding electromagnetic compatibility (see section 9).

Hazard to persons due to improper handling

- Use the device only for its intended purpose.
 - Never use the device for aspiration in the low vacuum range (e.g., thoracic drainage).
 - If aspiration procedures are carried out too often, slight bleeding and pain may occur.
 - When using the power supply unit, ensure that it is first connected to the aspirator and then to the mains supply (100 V to 240 V AC).
 - Disconnect the power supply unit from the mains supply in exactly the opposite sequence (first remove the power supply unit from the mains supply (100 V to 240 V AC) and then remove it from the aspirator).
 - Never touch parts of other devices than ME devices in the patient environment and the patient simultaneously.
-

Safety defects due to improper accessories and spare parts

The use of accessories and spare parts, detachable parts, or materials that are not recommended by Löwenstein Medical and specified in the instructions for use may compromise the safety and function of the device. Any warranty is excluded for damage caused by using accessories and spare parts not recommended or by improper use.

- Use only original accessories and spare parts. (see section 10)
-

Warning of safety defects due to improper connections to the ME system

The connection of the ME system with other devices, fixtures, or pieces of equipment not recommended by Löwenstein Medical and not specified in the instructions for use may compromise the safety and function of the ME system. Any warranty is excluded for damage caused by connecting devices, fixtures, or pieces of equipment not recommended with the **secretion aspirator Allegra M30** or by improper use.

- Connect only recommended original parts with the **secretion aspirator Allegra M30**.

Damage to the device due to heat

- Do not cover the power supply unit.
- Keep the aspirator, the mains cable, and the power supply unit away from other heat sources.

Damage to the device due to improper handling

- Never aspirate flammable, corrosive or explosive liquids or gases.
- Do not drop the device.
- Prior to each use, check the housing for any damage and do not operate the device if damage is apparent.
- Prior to each use, check all components that are subject to wear and damage to ensure that they are in perfect condition and proper operation of the device is guaranteed. If this is not the case, replace the parts immediately.

Possible physiological effects and unobvious risks

- Select the suction setting as appropriate for the respective patient and medical indication to avoid personal injuries.
- Always place the device upright on a sturdy and flat, incline-free base. Make sure that the device cannot be knocked over or fall and strike nearby persons.
- Never use the **secretion aspirator Allegra M30** to aspirate solids or inhaled foreign objects out of the throat. This could lead to the object being pushed deeper into the airways potentially causing a total obstruction.
- Operate the device only with the inserted vacuum generator and closed lid. The lid above the vacuum generator may be opened only by adequately trained service staff! Likewise, the vacuum generator may be replaced only by adequately trained service staff.
- Other devices, examinations, or treatments may be affected by the device. For this reason, always pay special attention to other devices and parallel examinations and treatments to detect any affect from the device as soon as possible.
- When using the device, always ensure adequate lighting so that all labels can be identified clearly.
- Small, detached parts could be inhaled or ingested. Therefore ensure that no unauthorized persons, children, or pets are near the device.
- Although the materials used have been tested for compatibility, in exceptional cases allergic reactions to the exposed materials on the device and its accessories may occur. This applies especially for contact injuries after prolonged exposure. In such cases, consult with a doctor immediately.

Known, identifiable or foreseeable conditions for medical care within a domestic environment

- Children and pets should be kept away from the device. Make sure that the device cannot be knocked over or fall and strike nearby persons.
- Prior to connecting the power supply unit, ensure that the device voltage corresponds to the domestic power supply.
- Ensure that lint and dust are immediately removed from the device and accessories in order to ensure unimpaired functionality. Furthermore, do not allow pests near the device because they may get inside and cause damage.

- Never expose the device and accessories to direct sunlight because the strong heat may cause functional impairments.
 - Some devices and sources that are commonly used in the domestic environment may cause potential faults in the device and accessories, e.g., fireplaces or radiant heaters (strong heat), and inhalers or steam kettles (excessive humidity). Do not operate such devices and sources in the vicinity of the **secretion aspirator Allegra M30**.
-

Damage to the device by ingress of liquids

The **secretion aspirator Allegra M30** has the IP classification IP22 against the ingress of liquids. Nevertheless, protect the device from moisture.

- Do not use the device near splashing water.
 - Do not use the device in damp rooms or while bathing or showering.
 - Never submerge the device in water or other liquids (including while it is switched off).
-

Contact of the interior of the device with liquids and solids

If the interior of the device comes into direct contact with liquids or solids, the device must be inspected by Löwenstein Medical or by an authorized service partner.

2.7 User requirements

The **secretion aspirator Allegra M30** may be operated and used only by instructed and appropriately trained persons. Familiarize yourself with the operation of the **secretion aspirator Allegra M30** prior to startup.

Training on how to handle the **secretion aspirator Allegra M30** is provided by Löwenstein Medical or an authorized distribution partner. Product training includes, among other things, an explanation of the design and function of the device, handling of the device, cleaning and optional disinfection, as well as the procedure to be followed for each new patient and for disposal.

2.8 Information on product liability

The liability for the operation of the device is transferred to the operator in the following cases:

- the device is used beyond the scope of its intended use;
- the device is not used in accordance with the instructions for use;
- the device is opened by unauthorized personnel;
- the safety seal is removed / broken;
- installation, settings, enhancements, processing or repairs are performed by unauthorized personnel;
- original accessories and spare parts have not been used.


The assembly of ME systems and modifications during their expected service life must be evaluated again based on the requirements of the applicable standards.

3 Product description

3.1 Complete view of the ME system


3.1.1 Secretion aspirator Allegra M30 with reusable secretion canister system

- A Reusable secretion canister system
- B Control panel
- C On / Off button
- D Socket for power supply unit (12 V DC)
- E Power supply unit
- F Mains cable
- G Suction tube
- H Carrying handle

Fig. 1

3.1.2 Secretion aspirator Allegra M30 with disposable secretion canister system

- A Disposable secretion canister system
- B Control panel
- C On / Off button
- D Socket for power supply unit (12 V DC)
- E Power supply unit
- F Mains cable
- G Suction tube
- H Carrying handle
- I Holder for external canister „Bag“

Fig. 2

en

3.2 Product contents

- The **secretion aspirator Allegra M30**
- These instructions for use
- Power supply unit incl. mains cable
- Charging instructions, multilingual
- "Used Medical Device" label and decontamination certificate
- Test report in accordance with IEC 62353
- Secretion canister system:
 - reusable secretion canister system (MW)** *OR* **disposable secretion canister system (EW)**
 - reusable secretion canister (1 l)
 - lid reusable secretion canister
 - suction tube
 - rinsing bottle (250 ml)
 - tube mount
 - disposable liner "OneWay"
 - external canister "Bag"
 - holder external canister "Bag"
 - suction tube
 - rinsing bottle (250 ml)
- Optional accessories (depending on the order, see section 10)

3.3 Product properties

The **secretion aspirator Allegra M30** is an especially handy, battery-powered medical aspiration device for mobile and stationary use.

The **secretion aspirator Allegra M30** has three suction settings (Low, Medium, High) and has a max flow rate of max. 30 l/min ± 2 l/min (see section 8 "Technical data").

The **secretion aspirator Allegra M30** is operated via the internal battery or via the supplied power supply unit, which can also be used to recharge the battery.

After 5 minutes of operation, the device switches off automatically. It is recommended to leave the device switched off for 5 minutes. However, the device can be switched on again immediately, as needed. Moreover, overheating protection prevents the battery from overheating by interrupting the charging process if the temperature of the battery is too high (e.g., due to unfavorable environmental conditions).

After it is switched on, the vacuum generator generates a vacuum, which is used to aspirate secretion (via a sterile aspiration catheter). The aspirate is directed away from the patient and collected in the secretion canister. If the secretion canister is full, the aspiration is stopped by means of overflow protection integrated into the lid of the secretion canister.

The desired suction setting can be selected by pressing the buttons on the control panel (Fig. 1(B) and Fig. 2(B), see section 2.2.2). The vacuum connection is established between the aspirator and secretion canister.

The device is charged or operated via the supplied 12 V power supply unit (Fig. 1(E) and Fig. 2(E)).

3.4 Information on the bacterial filter and the adapter for disposable secretion canister system

The **bacterial filter** or the **adapter for disposable secretion canister system** is located between the vacuum generator and the disposable secretion canister system.

The **bacterial filter** must be used when the **reusable secretion canister system** is connected. It is made of a hydrophilic (self-sealing) material and is effective against bacteria and viruses. It effectively protects the vacuum generator from overflow and contamination.

The **bacterial filter** is only intended for **reuse on one patient**.

The **bacterial filter** must be replaced each time a new patient uses the device (see section 4.2.2), in case of troubles (see section 6) or during **maintenance / repair**.

The **adapter for disposable secretion canister system** can be identified by the inscription "no filter / kein Filter" (Fig. 3) and may only be used with the disposable secretion canister system, as the adapter itself does not have a filter. The bacteria filter is located in the lid of the disposable liner "OneWay". This effectively protects the vacuum generator from overflow.

The **adapter for disposable secretion canister system** must be replaced each time a new patient uses the device (see section 4.2.2) or during **maintenance / repair**.


Fig. 3


Please observe the correct combination:

- Bacterial filter → reusable secretion canister system
- Adapter → disposable secretion canister system

3.5 Reusable secretion canister system for secretion aspirator Allegra M30

The reusable secretion canister system consists of the secretion canister, the lid for the secretion canister incl. overflow protection (see Fig. 3) and the suction tube with fingertip.

The overflow protection consists of a flexible plastic holder and a ball which is inserted into the plastic holder. The ball floats upwards in the holder when the canister is full and closes the transition from the upper part of the lid to the secretion aspirator. Thus, the suction process is interrupted. The buttons Low, Medium and High of the control panel light up red. Overflow of the **secretion aspirator Allegra M30** is thus prevented. The reusable secretion canister must be emptied and cleaned according to section 5.2.2.

- A Lid reusable secretion canister
- B Overflow protection with ball

Fig. 4

The reusable secretion canister system is for reuse on one patient.

3.6 Disposable secretion canister system for secretion aspirator Allegra M30

The disposable secretion canister system consists of the external canister "Bag", the holder for the external canister "Bag", the disposable liner "OneWay" and the suction tube. The disposable liner "OneWay" has an internal bacterial filter, carbon filter, and solidifier. The self-sealing bacterial filter integrated in the disposable liner "OneWay" is effective against microorganisms. It prevents overflow in the event of an operational error. If the liquid reaches this filter, aspiration is no longer possible and the vacuum setting buttons on the control panel illuminate continuously red. The aspiration process is discontinued. The disposable liner "OneWay" must be replaced.

The activated carbon filter in the disposable liner "OneWay" reduces the spread of odor.

Solidifier:

Disposable liners "OneWay" that are filled with aspirate can be transported and disposed of in a leak-proof manner with the aid of the solidifier. The aspirate solidifies after an average gelling time of 2 to 5 minutes (depending on the consistency of the secretion), irrespective of the aspiration intervals.

The **disposable liner "OneWay" and the suction tube** are for reuse on one patient.

Replace the disposable liner "OneWay" and suction tube according to the instruction in section 4.3.3 and 5.2.5.

3.7 Information on the battery

The charge state of the battery is indicated by the corresponding LED status indicator (for explanation, see section 4.3.1).

Prior to the initial startup of the **secretion aspirator Allegra M30**, it is urgently recommended to charge the battery completely and to repeat this after the initial uses.

The **secretion aspirator Allegra M30** is equipped with a lithium ion battery. Frequent short charging should be avoided. The typical service life of the battery is approximately 500 load cycles; thereafter approximately 70% of the primary capacity is still available.

The battery of the **secretion aspirator Allegra M30** is protected against total discharge; nonetheless, the above charging instructions should be observed. The battery is also protected against overheating during charging. If the battery temperature is exceeded during charging owing to improper ambient conditions, charging will be automatically interrupted temporarily to allow cooling. This feature ensures the safety of the battery and protects it from wear.

en

3.8 Warranty

Modifications of the ME device are not permitted.

The **secretion aspirator Allegra M30** is covered by warranty for 2 years. No extension or renewal of the warranty is granted in the event that warranty work is performed. The battery is covered by warranty for 6 months. Wearing parts are excluded from the warranty.

Löwenstein Medical is responsible for effects on safety, reliability and specific performance only if:

- original accessories and replacement parts are used,
- processing and repair are performed by authorized professionals or by Löwenstein Medical itself,
- the **secretion aspirator Allegra M30** is used and operated in accordance with the instructions for use and within its intended use.

All warranty claims shall be void if:

- the device is opened by unauthorized personnel,
- the safety seal is removed / broken,
- repairs are performed by unauthorized personnel,
- modifications are made to the device,

because the basic safety of the device can no longer be guaranteed in such cases and functional limitations may occur.

4 Operation

Before using the **secretion aspirator Allegra M30**, familiarise yourself with the basic safety instructions in section 2.6.

4.1 Setup and startup

The following section describes the operating elements, connections, and startup of the **secretion aspirator Allegra M30**:

4.1.1 Connecting the secretion aspirator Allegra M30

Check the power supply unit and mains cable for potential damage before each use and replace it if there is any damage.

Use the socket for the 12 V power supply unit (Fig. 1(D) and Fig. 2(D)) to connect the device to the mains power supply with the aid of the supplied power supply unit to charge or operate it as needed. Make sure that the device is positioned in such a way that it is possible to easily disconnect it later. First connect the power supply unit to the power supply socket of the **secretion aspirator Allegra M30** and then to the mains power supply.

Information on the permissible environmental conditions during operation can be found in section 8 “Technical data”.


Fig. 5 Connection of the **secretion aspirator Allegra M30** to patient and accessories

4.2 Startup

- Take the device and the accessories out of the packaging.
- Always place the device on a sturdy and flat, incline-free surface.
- Prior to the initial startup, it is essential to follow the safety instructions in section 2.6.
- Operate the device only with the inserted bacterial filter and connected secretion canister or with the disposable secretion canister and the associated adapter.
- Fully charge the battery prior to initial startup.
- Always have one backup bacterial filter (when used with the reusable secretion canister system) or an additional disposable liner "OneWay" (when used with the disposable secretion canister system) at hand, since it is absolutely essential for safe operation!
- Prior to each application, check that the reusable secretion canister is completely clean to avoid frothing.

en

4.2.1 Connecting the reusable secretion canister

3. Fold the moving part of the secretion canister lid (visioClick®) slightly upwards and place the secretion canister lid on the secretion canister.


Fig. 6


Fig. 7

4. Fold the moving part of the secretion canister lid (visioClick®) downwards so that the lateral hooks snap into place and the secretion canister is securely closed.

Fig. 8

5. Insert the secretion canister into the canister fitting provided on the device by guiding it over the rack on the bottom and pivoting it into the device until it snaps into place.


Fig. 9


Fig. 10

6. Connect the suction tube to the port in the middle of the secretion canister lid. Ensure that the tube is tight.


Fig. 11

4.2.2 Change of the bacterial filter or the adapter for the disposable secretion canister system

7. Remove the bacterial filter or the adapter for the disposable secretion canister system by moving the lever that is in the horizontal position counterclockwise a quarter turn so that the lever points vertically downwards.


Fig. 12


Fig. 13

8. When the lever is in this position, pull the bacterial filter or the adapter for the disposable secretion canister system out of the fitting.


Fig. 14

9. Insert the new bacterial filter or adapter for the disposable secretion canister system. For doing this, carry out steps 1 and 2 in reverse order.

4.2.3 Connecting the disposable secretion canister system


Malfunction due to collapsing disposable liner "OneWay"

A leak in the external canister „Bag“ or at the lid of the disposable liner „OneWay“ may cause air to flow into the external canister „Bag“. This may lead to the collapse of the disposable liner „OneWay“.

- Ensure that the lid of the disposable liner „OneWay“ is firmly connected to the external canister „Bag“.
- Ensure that all connections are firmly attached and properly connected.
- Ensure that the „Bag“ external canister is undamaged and the T-piece is firmly attached.

The original disposable secretion canister system consists of the external canister „Bag“, the holder for the external canister „Bag“, the disposable liner „OneWay“, the adapter for the disposable secretion canister system and the suction tube.

Connection designation

- A Vacuum connection
- B Patient connection

Fig. 15

10. Make sure that the adapter for the disposable secretion canister system is inserted correctly, otherwise the holder for the external canister „Bag“ cannot be inserted.

11. Insert the holder for the external canister „Bag“ into the canister fitting provided on the device by guiding it under the notch at the top and pushing it into the device via the rack at the bottom.


Fig. 16


Fig. 17

12. Remove the disposable liner „OneWay“ from the packaging and fully extend it.


Fig. 18

13. Place the disposable liner „OneWay“ in the reusable external canister „Bag“. Press the lid's edges firmly down to ensure proper sealing.


Fig. 19

14. Insert the bottom end of the T-piece of the external canister „Bag“ into the opening at the holder. Ensure that the bottom of the external canister „Bag“ is inserted into the canister fitting provided in the holder for the external canister „Bag“.


Fig. 20

15. Attach the prefitted connecting tube of the “OneWay” disposable liner to the upper end of the T-piece of the external canister „Bag“.


Fig. 21

16. Connect the patient connection of the disposable liner „OneWay“ (Fig. 15(B)) to the suction tube. Ensure that the tube is firmly attached.



Fig. 22

4.3 Operation

4.3.1 Control elements and indicators of the secretion aspirator Allegra M30

- A suction setting **High**
- B suction setting **Medium**
- C suction setting **Low**
- D LED status indicator for battery operation and battery charge level
- E LED status indicator for mains operation
- F On / Off button

en

Fig. 23

Symbol	Description	Meaning of display
	Button illuminates continuously green	Device is running, High vacuum setting is adjusted (-0.5 bar < vacuum ≤ -0.9 bar, flow rate 30 l/min)
	Button illuminates continuously red	Canister full / filter blocked (motor switches off automatically)
	Button illuminates continuously green	Device is running, Medium vacuum setting is adjusted (-0.2 bar < vacuum ≤ -0.5 bar, flow rate 25 l/min)
	Button illuminates continuously red	Canister full / filter blocked (motor switches off automatically)
	Button illuminates continuously green	Device is running, Low vacuum setting is adjusted (vacuum ≤ -0.2 bar, flow rate 20 l/min)
	Button illuminates continuously red	Canister full / filter blocked (motor switches off automatically)
	Display illuminates continuously green	Battery capacity 100% to 60%
	Display illuminates continuously orange	Battery capacity < 60% to 10%
	Display illuminates continuously red	Battery capacity < 10%, connect power supply unit immediately. Switch off the device after completion of the aspiration procedure for charging!
	Display blinks green	Battery is being charged.
	Display blinks red	Overtemperature protection, interruption of charging process until battery temperature returns to permissible range. OR Battery defective, if charging is not possible even if the device has cooled down. Operation only possible with connected power supply unit. Please contact your service partner!

Symbol	Description	Meaning of display
	Display illuminates continuously green	Power supply unit is connected, voltage available.
	Display blinks red	Incorrect power supply unit is connected, no charging of battery possible!
	Button illuminates continuously green	Device is switched on.
	Button illuminates continuously orange	Device still runs for 30 seconds before it switches off automatically. Afterwards the button will still be illuminated for 10 seconds.
	Button illuminates continuously red	Internal error, operation is not possible. Please contact your service partner!

4.3.2 Implementation of the aspiration

1. Switch on the device by pressing the On / Off button (Fig. 23(F)).
2. Select the desired suction setting by pressing the corresponding button (Low, Medium or High) on the control panel (Fig. 23(A, B, C)).
3. Perform the aspiration procedure as you were instructed by the specialised personnel. You can perform fine adjustments of the set flow rate via the additional air opening on the fingertip (Fig. 24).


Fig. 24

4. Switch off the device by holding the On / Off button (Fig. 23(F)) pressed for 1-2 seconds.

4.3.3 Replacement of the disposable liner „OneWay“

1. Switch off the **secretion aspirator Allegra M30** by holding the On / Off button (Fig. 2(C)) pressed for 1-2 seconds.
2. Separate the suction tube (Fig. 2(G)) from the patient connection (Fig. 15(B)).
3. Remove the external canister „Bag“ incl. the disposable liner „OneWay“ from the device.
4. Separate the pre-assembled connecting tube of the “OneWay” disposable liner at the upper end of the T-piece of the external canister „Bag“.

5. Remove the disposable liner „OneWay“ from the reusable external canister „Bag“.
6. Dispose of the disposable liner „OneWay“ in a proper manner (see section 7.3 “Disposal”).
7. Insert a new disposable liner „OneWay“ in the reusable external canister „Bag“ as specified in section 4.2.3. Make sure that the lid of the disposable liner is fitting tightly onto the external canister.
8. Attach the external canister “Bag” to the **secretion aspirator Allegra M30** according to section 4.2.3. Make sure that the lower T-piece is firmly seated in the holder for external canister “Bag”.
9. Plug the suction tube onto the patient connection (Fig. 15(B)) of the disposable liner „OneWay“.
10. Switch on the **secretion aspirator Allegra M30** by pressing the On / Off button (Fig. 2(C)).

4.4 Decommissioning

4.4.1 secretion aspirator Allegra M30 with reusable secretion canister

Observe the information on health risks when handling infectious or pathogenic germs in section 5.

5. Switch off the device after the aspiration by holding the On / Off button (Fig. 1(C)) pressed for 1-2 seconds.
6. Disconnect the power supply unit from the mains supply (100 V to 240 V AC) if attached, and then remove the device plug from the aspiration device.
7. Remove the suction tube (Fig. 1(G)) from the secretion canister.
8. Dispose of the suction tube properly.

9. Remove the reusable secretion canister from the device by pressing the unlocking mechanism on the secretion canister lid and then pivoting the secretion canister out of the device.


Fig. 25


Fig. 26

10. Fold the moving part of the secretion canister lid (visioClick®) upwards to release the lateral hooks at the lid and remove the lid from the secretion canister.


Fig. 27

11. Empty and dispose of the reusable secretion canister and the lid.
12. Remove the used bacterial filter according to section 4.2.2.
13. Clean the surface of the device according to section 5.2.1.
14. Store the device in the shipping carton until next operation.

4.4.2 secretion aspirator Allegra M30 with disposable secretion canister system

Note the hints for health risks due to handling of infectious or pathogenic germs in section 5.

1. Switch off the **secretion aspirator Allegra M30** after the aspiration by holding the On / Off button (Fig. 2(C)) pressed for 1-2 seconds.
2. Disconnect the power supply unit from the mains supply (100 V to 240 V AC) if attached, and then remove the device plug from the aspiration device.
3. Remove the suction tube (Fig. 2(G)) from the patient connection (Fig. 15(B)).
4. Dispose of the suction tube properly.
5. Remove the external canister "Bag" and the disposable liner „OneWay" as described in section 4.3.3 up to and including step 5 and dispose of these.
6. Remove the holder for the external canister „Bag" from the device by grasping it at the bottom in the middle and pulling it out of the canister fitting over the rack on the bottom (Fig. 28). Dispose of the holder for the external canister "Bag".


Fig. 28

7. Dispose of the adapter for the disposable secretion canister system according to section 4.3.3.
8. Clean the surface of the device according to section 5.2.1.
9. Store the device in the shipping carton until the next use.

5 Maintenance

5.1 General information

Health risks due to handling of infectious or pathogenic germs and

Infectious and pathogenic germs in the aspirate cause health risks.

- Wear suitable disposable gloves while doing so.
 - The reuse of bacterial filters for multiple patients is prohibited for hygiene and safety reasons!
 - For each new patient, processing by Löwenstein Medical or by an authorized service partner is strictly required!
 - After each aspiration procedure, all components that have come into contact with the aspirate must be cleaned or disposed of.
 - Disinfections can be done optional.
 - The aspirate and the parts contaminated with aspirate must be disposed of properly.
-

Damage of the device due to improper cleaning agents

- Do not use solvent-based cleaning agents for cleaning.
- To avoid damage or deface plastic parts of the components do not use solvent-based disinfectants.

Follow the instructions for use provided by the manufacturers of the disinfectants used, particularly concerning the compatibility with materials and surfaces as well as the concentration information and residence time.

Following disinfectants can be used: "Sekusept® aktiv" for immersion disinfection and "Incidin® Liquid" for wipe disinfection.

- Frequent cleaning and disinfection procedures may result in minor discolorations of the plastic components of the housing. However these do not affect the function of the device.
-

Health risks due to handling of disinfectants

- The use of appropriate protective clothing during disinfection is recommended.
 - Follow the disinfectant manufacturer's instructions.
-

Possible bodily injury due to electric shock

- Switch off the device prior to cleaning.
 - Disconnect the power supply unit by unplugging it from the mains supply. Subsequently disconnect the power supply unit from the socket for the power supply of the **secretion aspirator Allegra M30**.
-

Cleaning and optional disinfection of the device and the accessories do not affect the expected service life.

5.2 Processing

5.2.1 Surface of the device

Cleaning

The device can be wiped with a damp, lint-free cloth until no impurities are visible.


Manual Disinfection (optional)

A wipe disinfection can be done if internal hygiene requirements provide this.

Use disinfectants for wipe disinfection according to section 5.1 and observe the manufacturer's instructions.

en

5.2.2 Reusable secretion canister system

- A Reusable secretion canister lid
- B Overflow protection with ball

Fig. 29

Preparation before processing

1. Remove the reusable secretion canister system as described in section 4.4.1, step 5.
2. Separate the lid for reusable secretion canister from the reusable secretion canister as described in section 4.4.1
3. Empty the reusable secretion canister and dispose of the aspirate properly.
4. If you use an automated disinfection, please remove the overflow protection incl. ball (Fig. 29(B)). The ball can then be disinfected manually.

Cleaning

1. Rinse all components of the reusable secretion canister under running water under usage of a lint-free cloth.
2. Let all components of the reusable secretion canister dry.

Manual Disinfection (optional)

1. Immerse the components of the reusable secretion canister in the disinfectant solution. Use disinfectants for immersion disinfection according to section 5.1 and observe the manufacturer's instructions.
2. Rinse the components thoroughly and let them dry.

Automated thermal disinfection (optional)

The reusable secretion canister and the lid for reusable secretion canister (without overflow protection) can also be autoclaved for 20 minutes at 121 °C.

Reassembly

1. Fix the overflow protection in the lid for reusable secretion canister (see Fig. 29)
2. Put the ball in the overflow protection.
3. Connect the lid for reusable secretion canister with the reusable secretion canister and connect the reusable secretion canister system to the device (see section 4.2.1).

Replace the reusable secretion canister system if it is used and disinfected frequently, after 30 autoclaving processes at the latest, or if there are visible damages or functional restrictions, or if the patient changes.

5.2.3 External canister „Bag“

Preparation before processing

Disconnect the disposable liner “OneWay” and dispose of it as described in section 4.3.3, step 1-6.

Cleaning

1. Rinse the external canister “Bag” under running water under usage of a lint-free cloth.
2. Let the external canister “Bag” dry.

Manual Disinfection (optional)

1. Immerse external canister “Bag” in the disinfectant solution. Use disinfectants for immersion disinfection according to section 5.1 and observe the manufacturer’s instructions.
2. Rinse external canister “Bag” thoroughly and let it dry.

Automated thermal disinfection (optional)

The external canister “Bag” can also be autoclaved for 20 minutes at 121 °C.

Reassembly

1. Connect an new disposable liner “OneWay” to the external canister “Bag”.
2. Connect external canister “Bag” to the holder for external canister “Bag” (see section 4.2.3).

Replace replacing the external canister “Bag” if it is used and disinfected frequently, after 30 autoclaving processes at the latest, or if there are visible damages or functional restrictions, or if the patient changes.

5.2.4 Holder for the external canister "Bag"

Clean the surfaces of the holder for the external canister 'Bag' regularly and disinfect them daily.

Cleaning

The holder for external canister "Bag" can be wiped daily in course of the cleaning of the device in even if it is fitted, with a damp, lint-free cloth until no impurities are visible.

Manual Disinfection (optional)

A wipe disinfection can be done if internal hygiene requirements provide this. Use disinfectants for wipe disinfection according to section 5.1 and observe the manufacturer's instructions.

Replace the holder for external canister "Bag" if the patient changes.

5.2.5 Suction tube

Cleaning

1. Rinse the suction tube under running water.
2. Wipe up the suction tube with a lint-free cloth.
3. Let the suction tube dry.

Manual Disinfection (optional)

1. Immerse suction tube in the disinfectant solution. Use disinfectants for immersion disinfection according to section 5.1 and observe the manufacturer's instructions.
2. Rinse the suction tube thoroughly and let it dry.

Replace the suction tube if the patient changes.

If the suction tube is reused on one patient, replace the suction tube weekly at the latest.

In the event of heavy contamination with sticky or encrusted secretion, the suction tube must be replaced immediately.

5.2.6 Adapter for the disposable secretion canister system

Cleaning

1. Rinse the adapter for the disposable secretion system under running water.
2. Wipe up the adapter for the disposable secretion system with a lint-free cloth.
3. Let the adapter for the disposable secretion system dry.

Manual Disinfection (optional)

1. Immerse adapter for the disposable secretion system in the disinfectant solution. Use disinfectants for immersion disinfection according to section 5.1 and observe the manufacturer's instructions.
2. Rinse the adapter for the disposable secretion system thoroughly and let it dry.

Replace the adapter for the disposable secretion system if the patient changes.

5.3 Further use of the device on a new patient

The **secretion aspirator Allegra M30** is suitable for further use. However, prior to passing the device on to other patients or persons, it has to be processed professionally. For this purpose, please give the **secretion aspirator Allegra M30** to Löwenstein Medical or to qualified authorized personnel. In this regard, please observe the instructions in section 7.1!

6 Troubleshooting

Malfunction	Display	Probable cause	Remedy
Device does not start.	  do not illuminate	Battery is discharged	Connect the power supply unit for charging or operating.
	 blinks red	Battery is discharged and impermissible input voltage	Connect the correct power supply unit, check the power supply.
	 blinks red	Internal error	Please contact service!
Device does not charge.	 blinks red	Battery temperature too high for charging	No action required, automatic resumption of charging when battery has cooled down
	and  illuminates green	Battery is defective if charging is not possible when the device has cooled down.	Please contact service!
	 blinks red	Impermissible input voltage	Connect the correct power supply unit, check the power supply.
Insufficient suction	 illuminates green / orange	Leak in the suction line or on the canister	Check the lid of the disposable liner „OneWay“ or the reusable secretion canister lid and tube connections for tight fit and no leaks.
	and / or  illuminates green		
	 selected suction setting and  illuminate green	Bacterial filter is blocked	Replace the disposable liner “OneWay” or the bacterial filter

en

Malfunction	Display	Probable cause	Remedy
No suction	 continuously illuminate red	Ball of the overflow protection blocks the aspiration lid inlet (reusable secretion canister system)	Check fluid level in the reusable secretion canister, empty the reusable secretion canister.
		disposable liner „OneWay“ is full.	Replace the disposable liner „OneWay“.
	 selected suction setting and  illuminate green	Bacterial filter is blocked.	Replace the disposable liner „OneWay“ or the bacterial filter.
		Suction tube is blocked	Clean or change the suction tube
		Adapter for disposable secretion canister system or bacterial filter is not inserted	Inserting adapter for disposable secretion canister system or bacterial filter
Device switches off.	 blinks red	Internal error	Please contact service!

Contact Löwenstein Medical or your service partner if a malfunction cannot be corrected by the described measures.

7 Transport, storage, and disposal

7.1 Transport / Return

Löwenstein Medical offers its partners and customers fast and professional processing as well as required testing services.

Before the **secretion aspirator Allegra M30** is shipped to Löwenstein Medical, it must be cleaned and disinfected. Please follow the instructions in section 5.2! All disposable products must be disposed of correctly. The power supply unit and the mains cable must be enclosed with the return shipment for inspection. Please affix the supplied "Used Medical Device" label to the original shipping carton! Please notify Löwenstein Medical of the return in advance.

en

7.2 Storage

Store the **secretion aspirator Allegra M30** in accordance with the instructions in the technical data (section 8)!

Charge the battery of the **secretion aspirator Allegra M30** before storing the device. This ensures that it is operational at all times.

If the **secretion aspirator Allegra M30** is not used for a longer period of time (approximately 10 months), the battery must be fully recharged again!

7.3 Disposal

- The disposal of the device and the accessories must be carried out in a proper manner.
- Decontaminate the device and the accessories prior to disposal.
- According to EU Directives 2012/19/EU Waste Electrical and Electronic Equipment (WEEE) and 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS II), the device must not be disposed of in domestic waste.
- The device is registered in the national register for waste electric equipment (EAR) as small electrical device and can be turned into a collection station nearby for disposal.
- The device may also be disposed of via Löwenstein Medical or the service partner.
- Outside the EU: Follow the country-specific disposal requirements!

8 Technical data

Model name	Allegra M30
Flow rate* (measured at suction tube nozzle)	30 l/min ± 2 l/min (high flow)
Vacuum	max. -90 kPa (~ -900 mbar; -675 mmHg) (high vacuum) Conversion factor: 1 kPa ~ 7.5 mmHg; 1 kPa ~ 10 mbar
Canister	Reusable secretion canister (1 l) disposable liner „OneWay“ (1 l)
Suction tube (item number: 100713-2)	PVC, length 150 cm, vacuum stable up to -0.9 bar item number: 100713 inner diameter 6 mm, wall thickness 2 mm item number: 100765 inner diameter 7 mm, wall thickness 1.5 mm
Nominal voltage of the power supply unit	In: AC 100 – 240 V ~ 50-60 Hz / 1.5 A Out: DC 12 V / 5.0 A
Maximum load current	5.0 A
Permissible input voltage	12 V
Power consumption at 12 V	60 W
Degree of protection pursuant to IEC 60601-1	Type BF
Risk classification pursuant to 93/42/EEC, IX	IIa
Protection class pursuant to IEC 60601-1	II
IP degree of protection	IP22
CE marking	CE0494
Sound emission	max. 60 dB(A)
Ambient conditions	<i>Transport / Storage</i> Ambient temperature: -25 °C to +70 °C Relative humidity: up to 90%, non-condensing Air pressure: 700 hPa to 1060 hPa
	<i>Operation</i> Ambient temperature: +5 °C to +40 °C Relative humidity: 15% to 90%, non-condensing Air pressure: 700 hPa to 1060 hPa
Battery, rechargeable	at least 10.8 V; lithium ion battery
Charging time if battery is empty	at least 2 h
Energy of battery pack	<80 Wh
Power supply unit	GTM91099-6015-3.0-T2
Dimensions of device (H x W x D)	213 mm x 210 mm x 124 mm
Weight (device)	2.8 kg
Operating time	5 minutes ON; 5 minutes OFF
Runtime in battery operation	approximately 45 minutes, depending on the runtime of the motor
Expected service life	5 years
Item number (REF)	1510215h1 (with reusable secretion canister system) 1510220h1 (with disposable secretion canister system)

* The information provided may differ depending on the altitude above sea level, the prevailing air pressure, and the air temperature.

9 Electromagnetic compatibility

The **secretion aspirator Allegra M30** has been tested and found to comply with the limits specified in EN 60601-1-2:2015.

These limits and test levels are designed to provide reasonable assurance regarding electromagnetic disturbance when the device is used in its intended environment.

For more information on electromagnetic immunity and electromagnetic radiation, please contact Löwenstein Medical.

10 Order information

Item number	Description	PU
1510214hl4	Reusable secretion canister Allegra M30 (1 l)	1
1510214hl5	Lid for reusable secretion canister Allegra M30	1
1510214hl31	Reusable secretion canister system Allegra M30	1
1510214hl7	Universal holder Allegra M30	1
1510214hl18	Roller stand without device holder	1
1510214hl11	Device holder for roller stand Allegra M30	1
15102104hl8	Bag for Allegra M30	1
1510214hl9	Rinsing bottle (250 ml) Allegra M30	1
1510214hl10	Tube mount Allegra M30	5
1510214hl17	Car connection cable	1
100713	Suction tube with funnels and fingertip (inner diameter 6 mm, wall thickness 2 mm) (unsterile)	10
100765	Suction tube with funnels and fingertip (inner diameter 7 mm, wall thickness 1.5 mm) (unsterile)	10
1510214hl22	Holder rinsing canister (1 l) Allegra M30	1
15102104hl20	Set rinsing canister (1 l), lid, holder Allegra M30	1
1510214hl23	Lid for rinsing canister (1 l) Allegra M30	1
1510214hl22	Holder external canister „Bag“ (1 l) for Allegra M30	1
1510214hl21	Disposable secretion canister system for Allegra M30	1
1510214hl28	External canister „Bag“ (1 l)	1
100002	Disposable liner „OneWay“ (1 l)	60
1510214hl3e	Bacterial filter Allegra M30	10
1510214hl25	Adapter for disposable secretion canister system (1 l) Allegra M30	1

11 Publishing information

Created and published by:

Distributed by: Löwenstein Medical SE & Co. KG
Arzbacher Straße 80
56130 Bad Ems
GERMANY

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Manufacturer: ASSKEA GmbH