SOMNOvent CR

CS-Therapiegerät mit autoTriLevel-Prinzip/CS therapy device with autoTriLevel principle/Appareil de thérapie CS à principe autoTriLevel

SOMNOvent CR WM 24720

SOMNOvent CR mit/with/avec

SOMNOclick 300 WM 24785

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General information on setting therapy parameters

These instructions explain the setting of therapy parameters SOMNOvent CR directly on the device. Alternatively, the parameters can be set with the remote setting SOMNOadjust and the PC software WEINMANNsupport. These instructions are intended for medically trained experts only. Observe the operating instructions enclosed with each device!

How to access adjustment mode

1. The device must be plugged into an electrical socket and switched off.
2. Press the menu key , until "0" flashes in the display, and then press the humidifier key  as well. After about 4 seconds, the first value to be set will appear.
3. Increase the value to be set with the softstart button , or reduce it with the humidifier key . If the key is held down, the displayed time will change more rapidly. Press the menu key briefly  to call up the next value to be set.
4. To leave the setting mode, press the On/Off key/ESC . If no key is pressed within 3 minutes, the device will automatically leave the settingmode.

Values to be set:

Parameter		Value Range
Lower limit EEPAP	P_{\min} EEPAP	6 hPa to 17 hPa
Upper limit IPAP	P_{\max} IPAP	(Set lower limit EEPAP + 3 hPa) to 20 hPa
Upper limit EEPAP	P_{\max} EEPAP	EEPAP _{min} - IPAP _{max} - 3 hPa
Minimum respiratory frequency	$\min F$	6-10/min RuTo
Autostart	RuTo	on OFF
Mask test	P tEst	--, 8, 12, 16, 20 hPa
Softstart pressure	P 	4 hPa to lower limit EEPAP
Softstart time	\min 	5 to 30 minutes
Humidifier level		1 to 6
Time		h/min
Date		D/M/Y
Delete data	CLEAR dAtA	

Deleting therapy data

Therapy data is deleted e.g. on handing the device over to another patient. After the date, "CLEAR DATA" appears in the display. To delete the data, press the humidifier key (☹), until "CLEAR" flashes. Deletion cannot be undone.


General locking of Softstart

If the Softstart pressure is set to < 4 hPa, "--" appears in the display. Softstart is now locked and cannot be activated by the patient.

Connection to WEINMANNsupport and a polysomnograph

You can connect the SOMNOvent CR to WEINMANNsupport and the analog inputs of a PSG at the same time. For this, the Y-diverter WM 24972 is required, which is plugged into the front side socket of the SOMNOvent CR.

The connector cable to the PC for the WEINMANN*support* and the PSG connector cable can also be connected to this.

For this, the following accessories are available:

Y Part	Order Number
PSG connector cable without plug	WM 23976
Connector cable for Weinmann PSG	WM 24037
Connector cable for Respirationics / H+L-PSG	WM 23904
Connector cable for Viasys-PSG (Jäger)	WM 23908
Connector cable for MAP-PSG	WM 23929
Connector cable for Stimotron-PSG	WM 23941
Y-diverter	WM 24972
WEINMANN <i>support</i>	WM 93305

Two of the three parameters pressure, flow and leakage can be read out of the device directly to a PSG. The selection is made in WEINMANN*support*. Because of this, a pneumotachograph does not need to be used for therapy monitoring. A list of authorized pneumotachographs can be obtained from Weinmann.

Course of therapy

SOMNOvent CR saves the following data in detail for up to 10 hours:

- Actual pressure
- Target pressure (for IPAP, EPAP, EEPAP)
- Flow
- Relative respiratory volume per minute
- Leakage
- Respiratory events (e.g. apneas, hypopneas, RERA, flow limitations, snoring)
- Status (inspiration, expiration, artefacts, leakage, mask open)

The data for pressure, leakage and respiratory events (apnoeas, hypopnoeas, flow restrictions, snoring etc.) for the previous 30 therapy nights are saved in a time-compressed form.

Readout and evaluation of this data is only possible with WEINMANN*support* readout software.

For further details of the procedure for the readout and evaluation of the course of therapy and the setting of the parameters, please refer to the operating instructions WEINMANN*support*.

Readout of therapy data

The procedure for readout of the therapy data from the device is described in the operating instructions WM 66321.

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Geräte für Medizin GmbH+Co.KG