**DE** Gebrauchsanweisung **EN** Instructions for Use **FR** Mode d'emploi  
**IT** Istruzioni d'uso **ES** Manual de instrucciones

## Silentflow 3

Ausatemsystem  
Exhalation system  
Système expiratoire  
Espiratore  
Sistema de espiración

**LÖWENSTEIN**  
medical

<b>Deutsch</b> . . . . .	<b>3</b>
<b>English</b> . . . . .	<b>18</b>
<b>Français</b> . . . . .	<b>33</b>
<b>Italiano</b> . . . . .	<b>48</b>
<b>Español</b> . . . . .	<b>63</b>

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# Contents

<b>1</b>	<b>Intended use</b>	<b>20</b>
<b>2</b>	<b>Safety</b>	<b>20</b>
2.1	Safety information	20
2.2	General information	22
2.3	Warnings in this document	22
<b>3</b>	<b>Operation</b>	<b>23</b>
3.1	Assembling the exhalation system	23
3.2	Connecting the exhalation system	24
3.3	Dismantling the exhalation system	25
<b>4</b>	<b>Hygiene treatment and disposal</b>	<b>25</b>
4.1	Cleaning the exhalation system	26
4.2	Disinfecting the exhalation system	27
4.3	Disposal	28
<b>5</b>	<b>Faults</b>	<b>28</b>
<b>6</b>	<b>Technical specifications</b>	<b>28</b>
<b>7</b>	<b>Appendix</b>	<b>29</b>
7.1	Markings on the packaging	29
7.2	Characteristic pressure/flow curve	31
7.3	Materials	31
7.4	Warranty	32

7.5 Declaration of conformity .....	32
-------------------------------------	----

# 1 Intended use

The Silentflow 3 exhalation system is an accessory that removes CO<sub>2</sub> from the nasal mask and from the circuit. The exhalation system Silentflow 3 is used in combination with a therapy device in the treatment of sleep apnea and for non-invasive and invasive ventilation of patients with respiratory insufficiency.

Use of an exhalation system is mandatory for non-invasive ventilation with nasal masks and for invasive ventilation via a leakage circuit system.

## 2 Safety

### 2.1 Safety information

#### **Risk of injury from re-inhaling CO<sub>2</sub>!**

If the exhalation system is used incorrectly, CO<sub>2</sub> may be re-inhaled.

- ⇒ Do not seal off any of the openings of the exhalation system.
- ⇒ Only use the exhalation system within the stated therapy pressure range.
- ⇒ Only use the exhalation system when the therapy device is running.

#### **Risk of injury from anesthetic gases!**

Anesthetic gas may escape through the exhalation system and put third parties at risk.

- ⇒ Never use the exhalation system during anesthesia.

**Risk of injury from incorrect use!**

Incorrect use of the exhalation system may jeopardize the success of therapy.

- ⇒ Note the contraindications for mask and therapy device.
- ⇒ Do not use vented masks.
- ⇒ Ensure sufficient CO<sub>2</sub> wash-out of the selected system if using a non-vented full face mask.
- ⇒ Connect the exhalation system directly to the mask and use in the closest possible proximity to the patient.
- ⇒ Use the exhalation system only for a single patient.
- ⇒ Ensure that the connections to the mask and to the circuit are firm and air-tight. The connections may not come loose as a result of therapy pressure.

**Risk of injury from lack of cleaning!**

The exhalation system may be contaminated and this can put the patient at risk.

- ⇒ Clean the exhalation system before using it for the first time (see ["4.1 Cleaning the exhalation system"](#), page 26).
- ⇒ Clean the exhalation system regularly.
- ⇒ When selecting the detergent, consider potential allergies.

**Risk of injury from parts breaking off!**

The exhalation system may have loose parts or parts which come off and these may place the patient/user at risk.

- ⇒ Perform a visual inspection before using for the first time and after every cleaning operation.

## 2.2 General information

- In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority
- Silentflow 3 is compatible with all masks and circuits with an ISO female 22 mm tapered connection.

## 2.3 Warnings in this document

Warnings indicate information relevant to safety.

Within procedures, you will find warnings in front of a step which presents a risk to persons or objects.

### 2.3.1 Warning

Indicates an unusually significant hazardous situation. If you ignore this instruction, severe irreversible or fatal injuries may result.


### 2.3.2 Info

Indicates useful information within procedures.


## 3 Operation

### 3.1 Assembling the exhalation system

1. Push circuit connection (3) into mask connection (1) until the connecting clips (2) engage.

## 3.2 Connecting the exhalation system


1. The marking on the exhalation system must face the patient.
2. Push mask connection (**1**) into the rotating sleeve of the mask.
3. Push circuit connection (**3**) into the circuit.

### 3.3 Dismantling the exhalation system

1. Disconnect the exhalation system from the mask and from the circuit.

2. To release the circuit connection from the mask connection, squeeze the connecting clips together.

## 4 Hygiene treatment and disposal

### **WARNING**

**Residual risk of infection from inadequate cleaning!**

Contamination may put the patient at risk.

- ⇒ For patients with a compromised immune system or particular background of illness, disinfect the exhalation system daily following consultation with the physician.

## 4.1 Cleaning the exhalation system

1. Dismantle the exhalation system (see "3.3 Dismantling the exhalation system", page 25).
2. Clean the individual parts in accordance with the table below:

<b>Action</b>	<b>Daily</b>	<b>Weekly</b>
Wash individual parts in hot water (min. 30 °C) and mild detergent (1 ml detergent to 1 l water) for 15 minutes.	<b>X</b>	
Clean individual parts thoroughly with a cloth or a soft brush when washing. Or: Put individual parts in the dishwasher (max. 65 °C).		<b>X</b>

3. Wash immersed individual parts with a soft brush for at least 3 minutes.
4. Pay attention to all creases and cavities.
5. Rinse the individual parts in clean water.
6. Allow the individual parts to air-dry.
7. Perform a visual inspection.
8. If necessary: Replace a damaged exhalation system.
9. Assemble the exhalation system (see "3.1 Assembling the exhalation system", page 23).

Discolorations do not impair the function of the exhalation system.

## 4.2 Disinfecting the exhalation system

1. Clean individual parts before disinfection (see "4.1 Cleaning the exhalation system", page 26).
2. Immerse mask parts in Sekusept® Aktiv solution at a concentration of 2% for 15 minutes.<sup>1</sup>
3. Wash immersed individual parts with a soft brush for at least 3 minutes.
4. Pay attention to all creases and cavities.
5. Ensure that there are no air bubbles on the individual parts.
6. Rinse individual parts after disinfection for at least 1 minute with water (at least drinking water quality, sterile distilled water or fully demineralized water).<sup>2</sup>
7. Allow the individual parts to air-dry.
8. Assemble the exhalation system (see "3.1 Assembling the exhalation system", page 23).

<sup>1</sup>A bactericidal, levurocidal, tuberculocidal, mycobactericidal, sporicidal, fungicidal (clean cond.) and virucidal effect is achieved at this concentration and time to take effect.

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## 4.3 Disposal

You can dispose of the exhalation system in household waste or in line with hospital regulations.

## 5 Faults

Fault	Cause	Remedy
Draft in the eye	Exhalation system is not correctly aligned.	Align the exhalation system so that the marking points towards the patient.
Therapy pressure is not reached	Exhalation system is not correctly connected to the mask and circuit.	Check the connections to the circuit and to the mask.

## 6 Technical specifications

Product class to MDR (EU) 2017/745	Ila
Dimensions (H x D)	61.3 mm x 36 mm
Weight	13.2 g
Dead space volume	16 ml
Therapy pressure range	4 hPa - 40 hPa
Circuit connection/mask connection: Tapered connection to EN ISO 5356-1 (geometrical)	∅ 22 mm (male)
Temperature range Operation Transport and storage	+5 °C to +40 °C -20 °C to +70 °C

Flow resistance at 50 l/min at 100 l/min	0.12 hPa 0.44 hPa
Quoted two-figure noise emission value to ISO 4871: - Sound pressure level - Sound power level - Uncertainty factor	24 dB(A) 32 dB(A) 3 dB(A)
Service life	5 years
Useful life in a clinical setting	up to 7 days <sup>3</sup>
Useful life in a domestic setting	up to 6 months <sup>3</sup>

<sup>3</sup>The materials of the exhalation system deteriorate if exposed to aggressive detergents, for example. In individual cases, it may be necessary to replace the exhalation system sooner.

The service life of the exhalation system varies as a function of daily usage time and therapy pressure.

## 7 Appendix


### 7.1 Markings on the packaging

Symbol	Description
<b>CE 0197</b>	CE symbol (confirms that the product conforms to the applicable European directives/regulations)
<b>REF</b>	Order number

Symbol	Description
	Note temperature limits
	Keep out of sunlight
	Lot number
	Manufacturer
	Follow instructions for use
	Unique device identifier (uniform product code for medical devices)
	Use by date
	Indicates the product is a medical device

## 7.2 Characteristic pressure/flow curve

The characteristic pressure/flow curve shows leakage flow as a function of therapy pressure.

## 7.3 Materials

The exhalation system consists of polycarbonate (PC) and polypropylene (PP). Only use the exhalation system after consulting your physician if you are allergic to any of the materials.

The individual parts contain no latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

## 7.4 Warranty

Löwenstein Medical Technology gives the customer a limited manufacturer warranty on a new original Löwenstein Medical Technology product and any spare part fitted by Löwenstein Medical Technology in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Silentflow 3 exhalation system	6 months

## 7.5 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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