DE Gebrauchsanweisung 3 <i>LM Flowsensor Neo disposable</i>	NO Bruksanvisning 39 <i>LM Flowsensor Neo disposable</i>	HR Upute za uporabu 75 <i>LM Flowsensor Neo disposable</i>
EN Instructions for use 7 <i>LM Flowsensor Neo disposable</i>	SV Bruksanvisning 43 <i>LM Flowsensor Neo disposable</i>	HU Használati utasítás 79 <i>LM Flowsensor Neo disposable</i>
FR Mode d'emploi 11 <i>LM Flowsensor Neo disposable</i>	TR Kullanma talimatı 47 <i>LM Flowsensor Neo disposable</i>	SR Uputstvo za upotrebu 83 <i>LM Flowsensor Neo disposable</i>
ES Instrucciones de uso 15 <i>LM Flowsensor Neo disposable</i>	RU Руководство по эксплуатации 51 <i>LM Flowsensor Neo disposable</i>	LT Naudojimo instrukcija 87 <i>LM Flowsensor Neo disposable</i>
PT Manual de instruções 19 <i>LM Flowsensor Neo disposable</i>	PL Instrukcja użytkowania 55 <i>LM Flowsensor Neo disposable</i>	LV Lietošanas instrukcija 91 <i>LM Flowsensor Neo disposable</i>
ZH 使用说明 23 <i>LM Flowsensor Neo disposable</i>	CZ Návod k použití 59 <i>LM Flowsensor Neo disposable</i>	MK Упатство за употреба 95 <i>LM Flowsensor Neo disposable</i>
IT Istruzioni per l'uso 27 <i>LM Flowsensor Neo disposable</i>	SK Návod na použitie 63 <i>LM Flowsensor Neo disposable</i>	UA Посібник з експлуатації 99 <i>LM Flowsensor Neo disposable</i>
NL Gebruiksaanwijzing 31 <i>LM Flowsensor Neo disposable</i>	SL Navodila za uporabo 67 <i>LM Flowsensor Neo disposable</i>	EL Εγχειρίδιο χρήστη 103 <i>LM Flowsensor Neo disposable</i>
DA Brugervejledning 35 <i>LM Flowsensor Neo disposable</i>	RO Instrucțiuni de utilizare 71 <i>LM Flowsensor Neo disposable</i>	ID Panduan Pengguna 107 <i>LM Flowsensor Neo disposable</i>

REF 0217115v01

CE 0197

⚠️ WARNUNG

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Hersteller

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⚠ WARNING

Repreparation of disposable products can have a negative effect on the product properties. A change in the surface structure, for example, can lead to cracking or a reduction in the tear strength as well as to microbiological aggregation of, e.g., spores, allergens or pyrogens. Changing the material properties can lead to release of particles.

Danger of harm to the patient!

- Only use disposable products once and dispose of them properly afterwards.

⚠ WARNING

Defective flow sensors may not be used.

Danger of harm to the patient!

- Check the flow sensor for damage to the housing and the sensor wires before use.
- Dispose of the flow sensor in case of visible damage.

⚠ WARNING

Modifications to the product.

Danger of harm to the patient!

- All modifications to the flow sensor are prohibited.

NOTICE

- See instructions for use of the ventilator.
- The responsible organization is liable for ensuring the compatibility of the ventilator with all other parts that are used for connecting to the patients before use.
- The flow sensor must be calibrated prior to patient use.
- Change the flow sensor according to your hospital's infection control procedures or in light of secretion deposits or atomization of medication.
- Follow internal hospital procedures for infection control.

NOTICE

- Löwenstein medical provides no guarantee for the proper functioning of the disposable flow sensor if it is re-used.
- Every incident with the ventilator which leads to severe injury or death of the patient or to a potential public health hazard must be reported to the manufacturer and to the authorities responsible.

Intended use

The LM Flowsensor Neo disposable is a sensor for flow measurement. It is used in combination with a ventilator for ventilating of preterms and neonates as well as infants and children.

Restriction of the place of application

The flow sensor may not be used in areas with a strong magnetic field.

Permissible operating personnel

The product is operated by a doctor or by an individual trained and qualified specifically for this activity on the doctor's instructions, whereby every user must be trained on the product and familiar with the instructions for use and operation of the product.


Important for nebulization

The flow sensor must be removed prior to nebulization.

Scope of delivery

The flow sensor is delivered clean and ready for clinical use.

Connecting the flow sensor

- Connect the flow sensor (2) with the Y-piece (1) of the hose system and connect the patient-side ventilation hose (3) to the output of the flow sensor.
- Connect the connecting lead (4) of the flow sensor to the ventilator (please follow the detailed instructions in the instructions for use of the ventilator).

Calibrating the flow sensor

The flow sensor must be calibrated prior to patient use.

- Please follow the detailed instructions in the instructions for use of the ventilator to calibrate the flow sensor.
- Calibrate the flow sensor regularly, after switching on the ventilator and after each sensor replacement.

Technical data

Power supply (via flow sensor cable): 5 V DC

Measuring range: 0.1 l/min – 32 l/min

Accuracy: $\pm 8\%$

Operating conditions: Temperature 10°C – 40°C; relative humidity 15% – 95%.

Storage/transport conditions: Temperature -20°C – 55°C; relative humidity 20% – 75%, non condensing.

Spare parts

The following connecting leads are available for connecting the flow sensor to the respective ventilator:

- 0217012v01 Cable LM Flowsensor Neo, screwable
- 0739001v01 Cable LM Flowsensor Neo, push-pull

Disposal

Used flow sensors can be contaminated. Dispose of the product according to the regulations and infection control procedures of your clinic as well as the valid local waste disposal regulations. Dispose of the flow sensor cable in an environmentally friendly way as electrical scrap.

Symbols

Application part of type B (IEC 60417-5840)

Meets the requirements of the (EU) directive 2017/45 for medical products

See the instructions for use

Reference number

Do not use if the packaging is damaged

Batch no. XXXXX (XXXXX = consecutive number)

Do not re-use

Manufacturer

Use by

Representative in Switzerland

Do not dispose of with the domestic trash

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© 2024 Löwenstein Medical SE & Co. KG | Printed in Germany | All rights reserved
Löwenstein Medical Schweiz AG