DE Gebrauchsanweisung 3 <i>LM Flowsensor Neo reusable</i>	NO Bruksanvisning 48 <i>LM Flowsensor Neo reusable</i>	HR Upute za uporabu 93 <i>LM Flowsensor Neo reusable</i>
EN Instructions for use 8 <i>LM Flowsensor Neo reusable</i>	SV Bruksanvisning 53 <i>LM Flowsensor Neo reusable</i>	HU Használati utasítás 98 <i>LM Flowsensor Neo reusable</i>
FR Mode d'emploi 13 <i>LM Flowsensor Neo reusable</i>	TR Kullanma talimatı 58 <i>LM Flowsensor Neo reusable</i>	SR Uputstvo za upotrebu 103 <i>LM Flowsensor Neo reusable</i>
ES Instrucciones de uso 18 <i>LM Flowsensor Neo reusable</i>	RU Руководство по эксплуатации 63 <i>LM Flowsensor Neo reusable</i>	LT Naudojimo instrukcija 108 <i>LM Flowsensor Neo reusable</i>
PT Manual de instruções 23 <i>LM Flowsensor Neo reusable</i>	PL Instrukcja użytkowania 68 <i>LM Flowsensor Neo reusable</i>	LV Lietošanas instrukcija 113 <i>LM Flowsensor Neo reusable</i>
ZH 使用说明 28 <i>LM Flowsensor Neo reusable</i>	CZ Návod k použití 73 <i>LM Flowsensor Neo reusable</i>	MK Упатство за употреба 118 <i>LM Flowsensor Neo reusable</i>
IT Istruzioni per l'uso 33 <i>LM Flowsensor Neo reusable</i>	SK Návod na použitie 78 <i>LM Flowsensor Neo reusable</i>	UA Посібник з експлуатації 123 <i>LM Flowsensor Neo reusable</i>
NL Gebruiksaanwijzing 38 <i>LM Flowsensor Neo reusable</i>	SL Navodila za uporabo 83 <i>LM Flowsensor Neo reusable</i>	EL Εγχειρίδιο χρήστη 128 <i>LM Flowsensor Neo reusable</i>
DA Brugervejledning 43 <i>LM Flowsensor Neo reusable</i>	RO Instrucțiuni de utilizare 88 <i>LM Flowsensor Neo reusable</i>	ID Panduan Pengguna 133 <i>LM Flowsensor Neo reusable</i>

⚠️ WARNUNG

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Hersteller

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⚠ WARNING

The product is not delivered ready to use. The product may be contaminated after use on the patient.

Risk of infection!

- Prepare the product by validated procedures before using for the first time and after every use.
- Observe the national and institution-specific hygiene and preparation regulations.

⚠ WARNING

Signs of wear such as cracks, deformations, discolorations or peelings, for example, may occur on the prepared product.

Risk of infection!

- Observe the manufacturer specifications of the cleaning agent, disinfectant and preparation devices.
- Check for signs of wear after every product preparation and replace if necessary.
- Dispose of the flow sensor after using on the patient at the latest after 1 year or 52 preparation cycles.

⚠ WARNING

Modifications to the product.

Danger of harm to the patient!

- All modifications to the flow sensor are prohibited.

NOTICE

- See instructions for use of the ventilator.
- The responsible organization is liable for ensuring the compatibility of the ventilator with all other parts that are used for connecting to the patients before use.
- The flow sensor must be calibrated prior to patient use.
- Every incident with the ventilator which leads to severe injury or death of the patient or to a potential public health hazard must be reported to the manufacturer and to the authorities responsible.

Intended use

The LM Flowsensor Neo reusable is a sensor for flow measurement. It is used in combination with a ventilator for ventilating of preterms and neonates as well as infants and children.

Restriction of the place of application

The flow sensor may not be used in areas with a strong magnetic field.

Permissible operating personnel

The product is operated by a doctor or by an individual trained and qualified specifically for this activity on the doctor's instructions, whereby every user must be trained on the product and familiar with the instructions for use and operation of the product.

Important for nebulization

The flow sensor must be removed prior to nebulization.

Scope of delivery

The flow sensor is delivered individually packed and must be prepared before using for the first time and after every use.

Preparation (max. 52 cycles)

- Remove the flow sensor from the patient circuit and pull out the connecting lead if necessary.
- Rinse off the flow sensor immediately after use.
- Remove heavy contaminations from the sensor surface with a soft, moist cloth or a disposable wipe.

NOTICE

Observe the following restrictions for preparation of the flow sensor to avoid destroying the sensor wires:

- Never clean with compressed air or a hard water jet.
- Never clean in the ultrasonic bath.

Validated preparation procedure

Use only cleaning and disinfection devices (CDDs) which comply with ISO 15883 specifications for automatic preparation.

Automatic cleaning:

Pre-cleaning	with mains water	cold	for 2 min
Cleaning	with mains water and neodischer® MediZyme (Dr. Weigert) 0.5%	at 45°C	for 5 min
Rinse	with demineralized water	cold	
Disinfection	with demineralized water	at 90°C	for 5 min
Pre-drying	in the CDD	at 110°C	for 20 min

End drying:

End drying of the flow sensors is absolutely essential!

- Remove the flow sensor from the CDD and shake all remaining moisture out of the sensor in all directions (Caution! Do not knock the flow sensor).
- End dry the flow sensor in the drying cabinet.
Recommended drying time 30 min at 85°C

Sterilization (alternative to disinfection):

- Observe DIN EN ISO 17665 as a basic requirement.
- Sterilize the flow sensor with steam (as required) in an autoclave at e.g. $\geq 134^{\circ}\text{C}$ for 5 to 18 min.


Check:

- Carefully check that the surfaces of the flow sensor are clean and undamaged.
- Dispose of the sensor in case of visible damage.

Packing for keeping/storage:

- Pack the flow sensor individually in paper/film or in an equivalent packing of suitable size in accordance with DIN EN ISO 11607-1:2006 in connection with DIN EN 868.
- Do not use sterile goods packaging for disinfected products. (The packaging usually consists of mechanical protective packing and an outer packing, if necessary, regardless of the sterilization type.)

Connecting the flow sensor

- Connect the flow sensor (2) with the Y-piece (1) of the hose system and connect the patient-side ventilation hose (3) to the output of the flow sensor.
- Connect the connecting lead (4) of the flow sensor to the ventilator (please follow the detailed instructions in the instructions for use of the ventilator).

Calibrating the flow sensor

The flow sensor must be calibrated prior to patient use.

- Please follow the detailed instructions in the instructions for use of the ventilator to calibrate the flow sensor.
- Calibrate the flow sensor regularly, after switching on the ventilator and after each sensor replacement.

Storage and transport

- Store the flow sensor in a dry place protected from light.
- Observe the maximum storage times for sterile good in accordance with DIN 58953 Part 8.
- Avoid recontamination and damage during transport.

Technical data

Power supply (via flow sensor cable): 5 V DC

Measuring range: 0.1 l/min – 32 l/min

Accuracy: $\pm 8\%$

Operating conditions: Temperature 10°C – 40°C; relative humidity 15% – 95%.

Storage/transport conditions: Temperature -20°C – 60°C; relative humidity 10% – 100%, non condensing.

Spare parts

The following connecting leads are available for connecting the flow sensor to the respective ventilator:

- 0217012v01 Cable LM Flowsensor Neo, screwable
- 0739001v01 Cable LM Flowsensor Neo, push-pull

Disposal

Used flow sensors can be contaminated. Dispose of the product according to the regulations and infection control procedures of your clinic as well as the valid local waste disposal regulations. Dispose of the flow sensor cable in an environmentally friendly way as electrical scrap.

Symbols

Application part of type B (IEC 60417-5840)

Meets the requirements of the (EU) directive 2017/45 for medical products

See the instructions for use

Reference number

Do not use if the packaging is damaged

Batch no. XXXXX (XXXXX = consecutive number)

Sterilizable in the steam sterilizer (autoclave) at 134 °C

Manufacturer

Protect from light

Representative in Switzerland

Do not dispose of with the domestic trash

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© 2024 Löwenstein Medical SE & Co. KG | Printed in Germany | All rights reserved
Löwenstein Medical Schweiz AG