DE Gebrauchsanweisung für Module des Typs WM 100 MW
EN Instructions for Use for WM 100 MW Modules
FR Mode d'emploi pour modules du type WM 100 MW

prisma2CLOUD

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10.4	Warranty	45
10.5	Declaration of conformity	45

1 Introduction

1.1 Intended use

The WM 100 MW module is used to create a mobile radio connection between WM 100 TD, WM 110 TD, WM 120 TD therapy devices and the homecareONLINE web portal.

It is used to transmit compliance and therapy statistics from the therapy device to the homecareONLINE web portal.

The WM 100 MW module is intended for indoor domestic use only.

1.2 Function

The module is connected to the therapy device via the system interface. Due to the direct connection to the therapy device, the module does not have an on/off button.

The data communication is performed automatically.

Data communication is not possible during therapy.

For an explanation of the symbols in the display, please refer to the instructions for use of the therapy device.

2 Safety

Please read these instructions carefully. They form part of the devices described, and must be available at all times.

Use the unit for the designated purpose only (see "1.1 Intended use", page 26).

For your own safety and that of your patients, and in accordance with the requirements of Directive 93/42/EEC, please observe the following safety instructions.

2.1 Safety information

Warning

Risk of injury due to device or component malfunction!

A damaged device or damaged components may result in injury to the patient, user or bystanders.

- ⇒ Only operate the device and components if they are externally undamaged.
- ⇒ Only operate the device and components if the function check has been successfully completed.

Risk of injury if the device is operated outside the prescribed ambient conditions!

Use of the device outside the prescribed ambient conditions can result in failure to comply with tolerances, device failures, and injury to the patient.

- ⇒ Only operate the device within the prescribed ambient conditions (see "10.1 Technical data", page 43).
- ⇒ Only operate the device indoors.

Notice

Dirt in the device can cause material damage!

Dirt entering the device during transport can damage the device.

- ⇒ If the system interface port is not in use: Only transport the device if the port cover is closed.
- ⇒ Transport the device in the therapy device transport bag.

2.2 General information

- Use of third-party products may lead to functional failures and restricted fitness of purpose. Biocompatibility may also be compromised. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor original replacement parts recommended in the instructions for use are used.
- Repairs, servicing and maintenance should only be carried out by the manufacturer or by a technician expressly authorized by the manufacturer.
- The operator is responsible for ensuring the compatibility of the therapy device and all the connected components and accessories prior to the application with the patient. Only have modifications to the unit carried out by the manufacturer or by a technician expressly authorized by the manufacturer.
- Please observe the chapter on hygienic preparation in order to avoid infection or bacterial contamination (see chapter "Hygienic preparation").
- Also observe the respective instructions for use for the therapy device, the components, and the accessories.
- Please note that in the case of an existing LAN connection data communication via prismaCONNECT takes precedence, and no data are transferred via prisma2CLOUD.

2.3 Warnings in this document

Warnings are used to flag up safety-relevant information.

You will find a warning preceding any action that entails a hazard for persons or equipment.

Warnings consist of

- the warning symbol (pictogram),
- a signal word designating the hazard level,
- information about the hazard
- instructions for avoiding the hazard.

The warnings appear in three hazard levels depending on the degree of danger:

Warning!

Designates an extremely dangerous situation. Failure to observe this warning may lead to serious, irreversible, or fatal injury.

Caution!

Designates a dangerous situation. Failure to observe this warning may lead to minor or moderately serious injury.

Notice!


Indicates a hazardous situation. Failure to observe this warning may lead to damage to equipment.

Designates useful information relating to a particular action.

3 Product description

3.1 Overview


3-1 prisma2CLOUD module

No.	Designation	Description
1	System interface port cover	Covers the system interface port if no additional module is connected.
2	System interface connector	Used for connection to the therapy device or other modules.
3	Locking hooks	Affix the module to the therapy device or other modules.
4	Module unlocking button	Allows connection of the module to the therapy device and removal of the module from the therapy device.
5	Mounting holes	Accept locking hooks of other modules.

3.2 Labels and symbols

3.2.1 Labels on the module


3-2 Labels on the module

No.	Symbol	Description
Labels and symbols on the module		
1		Consult instructions for use
3		Mobile station

No.	Symbol	Description
Labels and symbols on the type plate		
2	TYPE: WM 100 MW	Type designation of prisma2CLOUD
	REF	Article number
	SN	Serial number
	SID	System identification number
		Manufacturer and date of manufacture
		Type of protection against electric shock: Protection class II device
		Do not dispose of device in household waste
	IP21	Degree of protection acc. to IP classification as per IEC 60529: Protected against ingress of solid foreign bodies ≥ 12.5 mm dia. Protected against vertically dripping water.
		Consult instructions for use
	CE mark (confirms that the product complies with the applicable European directives).	
	Device with non-ionizing electromagnetic radiation (mobile radio connection)	

3.2.2 Labels on the packaging

Symbol	Description
	CE mark (confirms that the product complies with the applicable European directives).
	Permissible storage temperature: -25°C to 70°C
	Permissible storage humidity: 15% to 93% relative humidity

Symbol	Description
	Protect the device against moisture.
	Fragile, handle with care
GTIN	Global trade item number
SN	Serial number
SID	System identification number

4 Preparation and operation


4.1 Connecting up the module

1. Place the therapy device on a level surface (e.g., a table).
2. Switch on the therapy device (see instructions for use of the therapy device).

You can also connect the module when the therapy device is switched off. Then switch on the therapy device.

3. Open the system interface cover on the right side of the therapy device.

4. Keep the unlocking button on the module depressed.
5. Slide the module's locking hooks into the mounting holes on the therapy device.
6. Push the module against the device until it clicks into place.


The  symbol appears on the therapy device's display.

Result Module is connected.

4.2 Removing the module

Requirement: There is no data communication underway.

1. Press the unlocking button on the module.

2. Remove the module sideways from the therapy device.
3. Close the system interface cover on the therapy device again.

Result The module is removed.
The  symbol disappears from the therapy device display.

5 Hygienic preparation

5.1 General information

- Wear suitable protective equipment for disinfection work.
- Please refer to the instructions for use supplied with the disinfectant used.

5.2 Intervals

Interval	Action
Weekly	Clean the modules
As necessary, e.g., following infectious diseases or severe soiling	Disinfect the modules

5.3 Hygienic preparation of the module

NOTICE

Damage to the device caused by ingress of liquids!

Ingress of liquids may damage the therapy device, components, and accessories.

⇒ Do not immerse the therapy device, components, or accessories in liquids.

Requirement The therapy device is switched off (see instructions for use for the therapy device).

1. Remove the module from the therapy device (see "4.2 Removing the module", page 36).
2. Carry out hygienic preparation of the module as specified in the following table:

Part	Cleaning	Disinfection	Sterilization
Housing	Wipe with a damp cloth using water or mild soap	Wipe-down disinfection (Recommendation: terralin® protect or perform advanced Alcohol EP)	Not permitted
Connection line			

3. Allow the module housing and in particular the electrical interfaces to dry well.
4. Connect the module up to the therapy device (see "[4.1 Connecting up the module](#)", page 34).
5. Switch on the therapy device (see instructions for use of the therapy device).

Result The module is hygienically prepared.

6 Function check

6.1 Intervals

Carry out a function check at regular intervals:

- After each hygienic preparation
- After each repair
- At least every 6 months

6.2 Performing a function check

1. Check whether the module is set up and connected correctly in accordance with the instructions for use.
2. Check the housing for cracks, damage, and heavy soiling.
3. Connect the module up to the device (see "[4.1 Connecting up the module](#)", page 34).
4. Check whether the  symbol is shown on the therapy device display.
5. If the module does not work properly or shows signs of damage: Contact your authorized dealer.

Result The function check is complete.

7 Faults

If you are not able to clear an error message with the aid of the table below, you should have the device repaired by your authorized dealer. To avoid serious damage, do not continue using the device.

Fault	Cause	Remedy
After connection, neither the green radio symbol  nor the gray radio symbol  appears on the therapy device display.	Contact problem or defective module.	Remove and then reconnect the module (see 4.2, p. 36), (see 4.1, p. 34). If the fault persists: Contact your authorized dealer.
	The module is not supported by the therapy device firmware.	Contact your authorized dealer and have the therapy device updated.
The gray radio symbol  appears on the therapy device display and does not turn green.	No communication possible between the module and therapy device.	Contact your authorized dealer.
	The module and the therapy device are not compatible.	

Please also observe the instructions for use of your therapy device.

8 Maintenance

The module is designed to have a useful service life of 6 years.

If the module is used as intended in accordance with these instructions for use, it does not require any maintenance.

If you identify faulty parts during the function check ([see "6 Function check", page 39](#)), contact your authorized dealer.


9 Storage and disposal

9.1 Storage

Store the device under the prescribed ambient conditions (see chapter "Technical data").

9.2 Disposal

9.2.1 Electronic waste

Do not dispose of the product in the household waste. Consult an authorized, certified electronic waste recycling company for proper disposal. You can find out their address from your environmental officer or from your local council.

The device packaging (cardboard box and inserts) can be disposed of as waste paper.

10 Appendix

10.1 Technical data

Specification	prisma2CLOUD
Product class according to 93/42/EEC	I
Dimensions W x H x D in cm	4.3 x 12.8 x 6.5
Weight in kg	0.170
Temperature range Operation	+5°C to +40°C
Storage	-25°C to +70°C
Permissible humidity during operation and storage	Rel. humidity 15% to 93%, non-condensing
Air pressure range	700 cmH ₂ O to 1100 cmH ₂ O (700 cmH ₂ O corresponds to normal pressure at a height of approx. 3000 m above sea level)
Power supply	12 V DC
Electrical output	1.5 VA*
Classification as per DIN EN 60601-11: Type of protection against elec. shock	Protection class II
IP classification as per ISO 60529: Degree of protection against harmful ingress of foreign bodies and water	IP21

Specification	prisma2CLOUD
Classification as per DIN EN 60601-1-11: Operating mode	Continuous operation
Electromagnetic compatibility (EMC) as per EN 60601-1-2 Radio interference suppression Radio interference immunity	Test parameters and limit values can be requested from the manufacturer if required. EN55011 B IEC 61000-4 Parts 2 to 6, Part 11, Part 8 IEC 61000-3 Parts 2 and 3
Radio module	 <p>The image shows a label for a Ublox LEON-G100 module. It includes the Ublox logo, a registered trademark symbol, and the text 'PP YY/WW' and 'xxS-xx'. Below this, it says 'Model: LEON-G100', '8XXXXX.XXXX.000', and 'IMEI: XXXXXXXXXXXXXXXXXX'. There is a QR code on the right. At the bottom, it shows 'CE 0682' and 'Made in Austria'.</p>
Radio standards	EN 301 511, EN 301 489-1, EN 301 489-7

*Data reached only in combination with the approved therapy device.

The right to make design modifications is reserved.

10.2 Scope of supply

prisma2CLOUD module

WM 29620

Part	Article number
prisma2CLOUD basic device	WM 29963
Instructions for use	WM 67950

10.3 Spare parts

Replacement parts can be ordered separately, if required. A current list of replacement parts is available on the manufacturer's website or from your authorized dealer.

10.4 Warranty

Starting from the date of purchase, Löwenstein Medical offers the customer a limited manufacturer's warranty on a new original Löwenstein Medical product or replacement parts installed by Löwenstein Medical in accordance with applicable warranty terms and conditions for the particular product and the warranty periods listed below. The warranty conditions can be downloaded from the manufacturer's website. We can also send you the warranty conditions on request.

In the event of a claim under warranty, please contact your specialist dealer.

Product	Period of guarantee
Löwenstein Medical devices, incl. accessories (excluding: masks), for sleep diagnostics, home ventilation, oxygen therapy, and emergency medicine	2 years
Masks, incl. accessories, batteries (unless otherwise stated in the technical documentation), sensors, hose systems	6 months
Disposable products	None

10.5 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG, Kronsaalsweg 40, 22525 Hamburg, Germany, declares herewith that the product complies fully with the respective regulations of the Medical Device Directive 93/42/EEC.

The complete text of the declaration of conformity is available on the manufacturer's website.

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