DE Gebrauchsanweisung für Module des Typs WM100MP EN Instructions for Use for WM100MP Modules FR Mode d'emploi pour modules tu type WM100MP

prismaPSG

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Contents

1	Introduction	28
1.1	Intended use	28
2	Safety	29
2.1	Safety information	29
2.2	General information	30
2.3	Warnings in this document	31
3	Product description	32
3.1	Overview of prismaPSG module	32
3.2	Markings and symbols	33
4	Preparation and operation	35
4.1	Connecting the prismaPSG module	35
4.2	Removing the prismaPSG module	36
4.3	Scaling on the PSG system	37
5	Hygiene treatment	38
5.1	General information	38
5.2	Intervals	38
5.3	Subjecting modules to a hygiene treatment	38
6	Function check	40
6.1	Intervals	40
6.2	Performing a function check	40
7	Faults	41
7.1	Faults when operating with prismaCONNECT	41
7.2	Faults when operating with prisma HUB	41
8	Servicing	42
9	Transport, storage, and disposal	43
9.1	Transport and storage	43
9.2	Disposal	43
10	Annex	44
10.1	Technical specifications	44

10.2	Scope of delivery	45
10.3	Accessories	45
10.4	Spare parts	45
10.5	Warranty	46
10.6	Declaration of Conformity	46

1 Introduction

1.1 Intended use

The WM100MP module is used to convert digital signals from the approved therapy devices to analog data. The analog signals present at the output of the WM100MP module can then be fed into a polysomnography system via the appropriate cabling.

The WM100MP module is used by health professionals in sleep labs for diagnosis and titration. The WM100MP module is designed solely for use indoors. The WM100MP module is not intended for use in a home environment.

2 Safety

Read these instructions for use carefully. They are a constituent part of the devices described and must be available at all times.

Use the device exclusively for the intended purpose described (see "1.1 Intended use", page 28).

For your own safety and the safety of your patients and in accordance with the requirements of (EU) Medical Devices Regulation 2017/745, follow the safety information below.

2.1 Safety information

- Warning* **Risk of injury from device or component malfunctions!**
A damaged device or damaged components may injure the patient, the user, and bystanders.
- ⇒ Only operate the device and its components if they are externally undamaged.
 - ⇒ Only operate the device and its components if a function check has been completed successfully.
- Risk of injury from operating the device outside the specified ambient conditions!**
Using the device outside the specified ambient conditions may lead to tolerances being exceeded and to device failure, injuring the patient.
- ⇒ Only operate the device within the specified ambient conditions (see "10.1 Technical specifications", page 44).
 - ⇒ Only operate the device indoors.
- Notice* **Material damage due to dirt in the device!**
Dirt which penetrates the device during transport may damage it.
- ⇒ Transport the device in the associated packaging.

2.2 General information

- If third-party articles are used, malfunctions and restricted fitness for use may result. Biocompatibility requirements may also fail to be met. In such cases, please be aware that any claim under warranty and liability will be void if neither the accessories nor the original spare parts recommended in the instructions for use are used.
- Have measures such as repairs, servicing, and maintenance work carried out by the manufacturer or by specialists expressly so authorized by the manufacturer.
- Connect only the therapy devices and modules approved in accordance with these instructions for use. Only medically-approved PSG systems must be connected to the prismaPSG module. Non-medical equipment should be positioned out of the patient's vicinity.
- The operator is responsible for ensuring the compatibility of the therapy device and of all the components or accessories connected to the patient before use.
Have modifications to the device carried out only by the manufacturer or an authorized specialist dealer.
- Follow the section on hygiene treatment ([see "5 Hygiene treatment", page 38](#)) to avoid infection or bacterial contamination.
- You should also follow the instructions for use of the therapy device, the components, and the accessories.
- In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

2.3 Warnings in this document

Warnings indicate information relevant to safety.

Within procedures, you will find warnings in front of a step which presents a risk to persons or objects.

Warnings consist of

- the warning symbol (pictogram),
- a signal word to indicate the level of hazard,
- information about the hazard and
- instructions on how to avoid the hazard.

There are three hazard levels for warnings depending on the degree of risk:

Danger!

Indicates an unusually significant hazardous situation. If you ignore this instruction, severe irreversible injuries or death will result.

Warning!

Indicates an unusually significant hazardous situation. If you ignore this instruction, severe irreversible or fatal injuries may result.

Caution!

Indicates a hazardous situation. If you ignore this instruction, mild or moderate injuries may result.

Notice!


Indicates a harmful situation. If you ignore this instruction, material damage may result.

Indicates useful information within procedures.

3 Product description

3.1 Overview of prismaPSG module



3-1 prismaPSG module

No.	Designation	Description
1	PSG connection	For connecting the PSG connecting cable to the prismaPSG module.
2	Analog output socket 8	Respiratory frequency
3	Analog output socket 7	Tidal volume
4	Analog output socket 6	Minute volume
5	Analog output socket 5	Actual pressure
6	Analog output socket 4	Setpoint pressure
7	Analog output socket 3	Obstruction level
8	Analog output socket 2	Unintentional leak
9	Analog output socket 1	Respiratory flow
10	PSG connecting cable	Connects the prismaPSG module to the prismaCONNECT module or prisma HUB.

3.2 Markings and symbols

The following symbols may be applied to the device, the device ID plate, accessories or packaging.

Symbol	Description
	Serial number
	Indicates the product is a medical device
	Unique device identifier
	Order number
	Manufacturer and date of manufacture
IP21	Degree of protection against contact with a finger. Product is protected against vertically-falling drips of water
	Electrical input
	Electrical output
	Do not dispose of the product in domestic waste
	Degree of protection against electric shock: Protection class II product
	Follow the instructions for use
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

Symbol	Description
	Permitted humidity range for transport and storage
	Permitted temperature range for transport and storage
	prismaCONNECT connection

4 Preparation and operation

4.1 Connecting the prismaPSG module


You can use the module with the following therapy devices:

- prisma SMART series (WM090TD) - via the prisma HUB module
- prismaLINE series (WM100TD) - via the prismaCONNECT module or prisma HUB
- prisma VENT series (WM110TD/WM120TD) - via the prismaCONNECT module or prisma HUB

The prismaPSG module (WM100MP type) is connected to the prismaCONNECT module (WM100MC type) or the prisma HUB module (WM090MC).

Requirement

- The prismaCONNECT module or prisma HUB is connected to the therapy device (see instructions for use of the prismaCONNECT module or prisma HUB).
1. Connect the plug of the PSG connecting cable supplied to the PSG connection of the prismaPSG module.

2. Connect the other plug of the PSG connecting cable to the module interface connection of the prismaCONNECT module or the prisma HUB.
3. Connect the prismaPSG module to a PSG. Preconfigured plugs and a universal cable are available for configuration ([see "10.3 Accessories", page 45](#)).
4. If necessary: Switch on therapy device.

Result prismaPSG has been connected.

4.2 Removing the prismaPSG module

1. Remove the plug of the PSG connecting cable from the module interface connection of the prismaCONNECT module or the prisma HUB.

Result The module has been removed.

4.3 Scaling on the PSG system

There are 8 analog sockets on the prismaPSG module. These are assigned fixed parameters:

Analog output socket	Signal	Lower phys. limit	Upper phys. limit
1	Respiratory flow	-250 l/min	250 l/min
2	Unintentional leak	0 l/min	100 l/min
3	Obstruction level	0 %	100 %
4	Setpoint pressure	0 hPa	40 hPa
5	Actual pressure	0 hPa	40 hPa
6	Minute volume	0 l/min	63 l/min
7	Tidal volume	0 ml	2000 ml
8	Respiratory frequency	0 bpm	63 bpm

To mark the assignment of the channels, cut out one of the paper strips on the last page of the instructions for use and push it into the control panel of the prismaPSG module.

5 Hygiene treatment

5.1 General information

- Wear appropriate safety gear for disinfection.
- Follow the instructions for use of the disinfectant used.

5.2 Intervals

Interval	Action
Weekly	Clean modules
If required, e.g., following infectious diseases or severe contamination	Disinfect modules

5.3 Subjecting modules to a hygiene treatment

NOTICE

Material damage due to ingress of liquids!

Ingress of liquids may damage therapy device, components, and accessories.

⇒ Do not immerse therapy device, components, and accessories in liquids.

1. Disconnect the prismaPSG module from the therapy device (see "[4.2 Removing the prismaPSG module](#)", page 36).
2. Subject the module to a hygiene treatment in accordance with the table below:

Part	Cleaning	Disinfection	Sterilization
Housing	Wipe down: Use water or mild detergent	Disinfect by wiping (Recommended: terralin® protect or perform advanced Alcohol EP)	Not permitted
Connecting cable			

3. Leave the module housing and in particular the electrical interfaces to dry thoroughly.
4. Connect the module to the device (see "4.1 Connecting the prismaPSG module", page 35).

Result The prismaPSG module has been subjected to hygiene treatment.

6 Function check

6.1 Intervals

Perform a function check at regular intervals:

- After every hygiene treatment
- After every repair
- At least every 6 months

6.2 Performing a function check

1. Check whether all the modules have been set up and connected properly in accordance with the instructions for use.
2. Check housing for cracks, damage, and severe soiling.
3. Connect the prismaPSG module ([see "4.1 Connecting the prismaPSG module", page 35](#)).
4. For operation with prismaCONNECT: Check whether the **PSG** symbol is shown in the display of the therapy device.
For operation with prisma HUB: Check whether the **prismaPSG** LED on the control panel of the prisma HUB is on.
5. Check whether a value matching what is shown in the display of the therapy device is displayed on the pressure channel being fed in.
6. If the module is not working properly or is damaged: Contact your specialist dealer.

Result The function check is complete.

7 Faults

If you are unable to remedy faults with the aid of the table, or in the event of unexpected operation or an incident, contact your specialist dealer to have the device repaired. To avoid exacerbating the damage, do not continue operating the device.

7.1 Faults when operating with prismaCONNECT

Fault	Cause	Remedy
The PSG symbol does not appear in the display of the therapy device	PSG connecting cable not plugged in or defective	Disconnect PSG connecting cable and plug it back in.
The PSG symbol appears in the display of the therapy device, but analog values are "0" or implausible.	Connection to PSG defective or non-existent.	Check plug connections to the PSG.

7.2 Faults when operating with prisma HUB

Follow the instructions for use of the prisma HUB module.

8 Servicing

The prismaPSG module is designed for a service life of 6 years.

If used in accordance with intended purpose in line with these instructions for use, none of the modules requires any maintenance.

If you find defective parts during the function check ([see "6 Function check", page 40](#)), contact your specialist dealer.

9 Transport, storage, and disposal

9.1 Transport and storage

Transport and store the device under the specified ambient conditions ([see "10.1 Technical specifications", page 44](#)).

9.2 Disposal

9.2.1 Electronic scrap

Do not dispose of the product in domestic waste. To dispose of it properly, contact a licensed, certified electronic waste disposal merchant. This address is available from your Environment Officer or from your local authority. The device packaging (cardboard and inserts) can be disposed of in paper recycling facilities.

10 Annex

10.1 Technical specifications

Specification	prismaPSG
Product class to MDR (EU) 2017/745	I
Dimensions W x H x D in cm	15 x 3 x 9.2
Weight in kg	0.2
Temperature range Operation	+5 °C to +40 °C
Transport and Storage	-25 °C to +70 °C
Permitted humidity for operation, transport, and storage	Relative humidity 15 % to 93 %, no condensation
Air pressure range	700 hPa to 1060 hPa, corresponds to an altitude of approx. 3000 m above mean sea level
Electrical rating	12 V DC
Power capacity	0.6 VA*
Outputs 1-8 Output voltage	0-1 V
Resolution	0.002 V
Socket assignment	1 channel per 3.5 mm jack bush, tip: signal, shaft: ground
Insulation against power supply	4 kV
Insulation against communication signals	2.5 kV
Insulation between the 8 channels	None

Specification	prismaPSG
Classification to IEC 60601-11: Type of protection against electric shock	Protection class II
Protection against damaging ingress of water and solids	IP21
Classification to IEC 60601-1: Operating mode	Continuous duty
Electromagnetic compatibility (EMC) to IEC 60601-1-2* Emissions Immunity	CISPR 11 Class A, IEC 61000-3 Parts 2 and 3 IEC 61000-4 Parts 2 to 6, 8, 11, 39

*Data apply only in combination with the approved device.

The right to make design modifications is reserved.

10.2 Scope of delivery

prismaPSG module

WM 29690

Part	Article number
prismaPSG basic device	WM 29691
PSG connecting cable to prismaCONNECT	WM 29696
Instructions for use for prismaPSG	WM 67960

10.3 Accessories

Part	Article number
PSG connecting cable, Heinen and Löwenstein (registered jack)	WM 35151
PSG connecting cable (registered jack to fit Transferbox and Transferbox 2, inputs 1-4)	WM 35152

Part	Article number
PSG connecting cable, 3.5 mm jack (e.g. Alice/Embla)	WM 35153
PSG connecting cable, 2.5 mm jack (e.g. Nihon Kohden)	WM 35154
PSG connecting cable, universal	WM 35155

10.4 Spare parts

You can order spare parts separately if necessary. A current list of spare parts can be obtained on the manufacturer's website or through your specialist dealer.

10.5 Warranty

Löwenstein Medical Technology gives the customer a limited manufacturer warranty on a new original Löwenstein Medical Technology product and on any spare part fitted by Löwenstein Medical Technology in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Devices including accessories (exception: Masks)	2 years
Masks including accessories, batteries (unless quoted differently in the technical documentation), sensors, circuits	6 months
Disposable products	None

10.6 Declaration of Conformity

Löwenstein Medical Technology GmbH + Co. KG, Kronsaaßweg 40, 22525 Hamburg, Germany, the manufacturer of the devices described in these instructions for use, hereby declares that the product complies with the relevant provisions of Medical Devices Regulation (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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