
LMT 68326SA0a 08/2024 AR, EN-US


CE 0197

**Manufacturer**  
Löwenstein Medical  
Technology GmbH + Co.KG  
Kronsaalsweg 40  
Hamburg, Germany 22525  
T: +49 40 54702-0  
F: +49 40 54702-461  
[www.loewensteinmedical.com](http://www.loewensteinmedical.com)







**CARA**  
Nasal Mask

LMT 68326SA0a

**LÖWENSTEIN**  
medical

**LÖWENSTEIN**  
medical

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with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor original spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

## 12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

### Contact information

Info@azeermed.com

+966551177970 Azeer Customer Support

+966126502921 Azeer Customer Support

www.azeermed.com

P.O. Box 132391 Jeddah 21382 KSA. Prince majed Street

## 1 Operation

The figures show the following steps for operating the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask

### For blind and partially-sighted users

An electronic version of the instructions for use is also available on the manufacturer's website.

## 2 Introduction

### 2.1 Intended use

The CARA mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

### 2.2 Contraindications

The mask must not be used for patients weighing < 30 kg.

The mask may not be used in the following situations: Necessity for immediate intubation, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: pressure points and acute injuries to the facial skin, skin allergies in the facial area, facial or nasopharyngeal deformities, acute pain in the facial area, cough reflex restricted or absent, claustrophobia, acute nausea.

Mask is not suitable for use with a nebulizer or as an anesthetic mask.

If you are not certain whether one of these situations applies to you, consult your healthcare professional. Please observe the contraindications in the instructions for use of your device.

### 2.3 Side effects

The following side effects may occur with use of the mask: Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

If these side effects occur, contact your healthcare professional.

### 2.4 Clinical benefit

Transfers the therapeutic efficacy of the therapy device to the patient.

## 3 Safety

**⚠ WARNING!** Unusually hazardous situation. If you ignore the following instructions, severe, irreversible or fatal injuries may result.

### Risk of injury due to damaged mask parts or those under strain!

- ⇒ Perform a visual inspection before every use and after every cleaning operation. Replace mask parts if necessary.
- ⇒ Note useful life (see section entitled "Technical data").

### Risk of injury from re-inhaling CO<sub>2</sub>!

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients unable to remove the mask themselves must be monitored by a nurse.
- ⇒ Do not close off exhalation systems.

### Risk of injury due to the use of oxygen!

Oxygen can accumulate in clothing, bed linen, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.
- ⇒ Observe the safety instructions in the instructions for use of the device.

### Risk of injury due to insufficient supply to the patient!

- ⇒ Activate low pressure/leakage alarms on device.
- ⇒ Use the correct mask size and check that it fits tightly.

### Risk of injury from inadequate cleaning!

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.
- ⇒ In case of a change of patient in the clinical environment: Follow the document *Mask reprocessing instructions* (see section entitled "Reprocessing").

**⚠ CAUTION!** Hazardous situation. If you ignore these instructions, mild or moderate injuries may result.

### Risk of injury due to allergies!

- ⇒ Be aware of the materials used in the mask (see section entitled "Materials"). Only use a mask following consultation with a healthcare professional.
- ⇒ When selecting the detergent, consider potential allergies.

## 4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	4	Rotating sleeve
2	Forehead support	5	Mask body
3	Elbow	6	Mask cushion

### Compatible devices

In some device combinations, the actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that the actual pressure in the mask corresponds to the therapy pressure. This adjustment should be made with the type of mask that will be used during therapy.

### Exhalation system

The mask has an integrated exhalation system. The exhaled air escapes through a gap.

## 5 Cleaning and reprocessing

### 5.1 Clean mask

1. Wash your hands before starting cleaning.
2. Dismantle mask (see Figure 4).
3. Clean the mask by hand (max. 30 °C, 1 ml mild detergent to 1 l water) in accordance with the following table:

Mask part	Frequency	Action
All mask parts	Daily	Soak and wash for 15 minutes and clean for 3 minutes using a soft cleaning brush.
Headgear	Weekly	Wash for 15 minutes.

**i** All parts can be washed in a dishwasher once a week (max. 70 °C, mild detergent, max. program length 90 minutes, top basket, separate rinse).

4. Rinse all parts with clean water.
5. Allow all parts to air-dry.
6. Carry out a visual inspection for cracks and deformations. Replace damaged parts. Discoloration is harmless.

7. Re-assemble mask (see Figure 5).

### 5.2 Reprocessing (clinical environment)

In the event of a change of patient, follow the instructions in the document entitled *Information on reprocessing*. The document can be found on the manufacturer's website. We will send you this document on request.

### 5.3 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

## 6 Troubleshooting

Fault	Cause	Action
Pain due to pressure on the face	Mask too tight.	Adjust headgear so that it is looser.
Draft in the eye	Full face mask fits too loosely.	Adjust headgear so that it is tighter.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask is not correctly adjusted.	Re-adjust mask.
	Mask cushion is damaged.	Replace mask cushion.
	Circuit is damaged.	Check circuit and the correct fit of circuit.

## 7 Technical data

Dimensions in mm (W x H x D)	
Size XS	72 x 112 x 70
Size S/M	72 x 117 x 70
Size M/L	72 x 118 x 70
Size XL	70 x 125 x 75
Weight	
Size XS	58 g
Size S/M	59 g
Size M/L	60 g
Size XL	69 g
Dead space volume	
Size XS	74 ml
Size S/M	76 ml
Size M/L	79 ml
Size XL	98 ml

Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (male)
Flow resistance at 50 l/min	0.2 hPa
at 100 l/min	0.9 hPa
Service life	5 years
Useful life	Up to 12 months <sup>1</sup>
Therapy pressure	4 hPa - 30 hPa
Quoted two-figure noise emission value according to ISO 4871:	
Sound pressure level	18 dB(A)
Sound power level	26 dB(A)
Uncertainty factor	3 dB(A)
Temperature:	
Operation	+5 °C to +40 °C
Transport and storage	-20 °C to +70 °C
Standards applied	EN ISO 17510: 2020
Product class to MDR (EU) 2017/745	IIa


<sup>1</sup> Useful life depends on cleaning and on the cleaning agents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. The use of the second mask cushion (optionally included) will not extend the useful life.

## 8 Materials

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

Headgear	CO (cotton), EL (elastane), P (polyester), PU (polyurethane)
Headgear clip	PA (polyamide)
Forehead support	PA (polyamide)
Mask body	PA (polyamide)
Mask cushion	SI (silicone)
Elbow	PA (polyamide)
Rotating sleeve	PA (polyamide)

## 9 Characteristic pressure/flow curve

## 10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Unique device identifier (unique product code for medical devices)
	Order number
	Indicates the product is a medical device
	Follow the instructions for use
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Lot number
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

## 11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance