
DE Gebrauchsanweisung **EN** Instructions for Use **FR** Mode d'emploi
NL Gebruiksaanwijzing **IT** Istruzioni d'uso **TR** Kullanma Kılavuzu
ES Manual de instrucciones **PT** Instruções de uso **RU** Инструкция по пользованию
DA Brugsanvisning **SV** Bruksanvisning **CS** Návod k použití **EL** Οδηγίες χρήσης

JOYCEclinic Full Face

Full Face Mask
vented, non vented

WM168110j 08/2023 DE, EN, FR, NL, TR, IT, ES, MX, PT, RU, DA, SV, CS, EL





CE 0197

LÖWENSTEIN
medical

**Löwenstein Medical
Technology GmbH + Co. KG**
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.de

LÖWENSTEIN
medical

12 Bestellinformation

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1 Operation

Use the template on the title page (see Figure 4) of these instructions for use to select the correct mask size. See the following figures for how to put on, adjust and remove the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask

2 Introduction

2.1 Intended use

The JOYCEclinic full-face mask is used as a connecting element between the patient and the therapy device for non-invasive positive pressure respiratory therapy in hospitals/clinics. It is for a single patient with respiratory insufficiency (> 30 kg) who is breathing spontaneously and who may be subjected to non-invasive ventilation with pressure support. The mask may only be used in combination with therapy devices which have the appropriate alarms and safety systems in case the device fails.

Use JOYCEclinic NV only in combination with therapy devices which have an active exhalation valve.

2.2 Contraindications

The mask may not be used, or may be used only with particular caution, if any of the following symptoms are present:

erosions and ulcerations, skin allergies, rashes on the face, pressure marks on the face, claustrophobia, anxiety, deformities of the face or nasopharynx, taking of medication which may induce vomiting, the necessity for immediate intubation, pronounced reflux, dysfunctional cough reflex. Please also observe the contraindications in the instructions for use for your therapy device.

2.3 Side effects

Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing. Masks may cause teeth, gums or jaw to ache.

3 Safety

Risk of injury from the supply of oxygen!

Oxygen can become deposited in clothing, bedlinen and hair. In conjunction with smoking, naked flames or electrical equipment, it can cause fires and explosions.

- ⇒ Do not smoke.
- ⇒ No naked flame.

Risk of injury from re-inhaling CO₂!

If the mask is used incorrectly, CO₂ may be re-inhaled.

- ⇒ Vented variant only: do not close off the exhalation system of the mask.
- ⇒ Only put on the mask for an extended period if the therapy device is running.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.

Risk of injury if the mask slips!

If the mask slips or falls off, the therapy is ineffective.

- ⇒ Monitor patients with restricted spontaneous respiration.
- ⇒ Activate low pressure/leakage alarms on the therapy device.

Risk of injury from anesthetic gases!

Aesthetic gas may escape through the exhalation valve and put third parties at risk.

- ⇒ Never use the mask during anesthesia.

Risk of injury from lack of cleaning!

The mask may show contamination which may put the patient at risk.

- ⇒ Clean the mask before using for the first time (see section entitled "Hygiene treatment").
- ⇒ Clean mask regularly.

4 Product description

4.1 Overview

A diagram of the individual parts can be found on the title page.

1. Headgear
2. Forehead cushion
3. Forehead support
4. Retaining ring
5. Elbow
6. Anti-asphyxia valve
7. Mask body
8. Mask cushion

9. Headgear clip
10. Quick-release cord
11. Rotating joint adapter (for JOYCEclinic FF vented only, not shown)

4.2 Compatible devices

With many combinations of device, actual pressure in the mask will not correspond to prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the combination of devices adjusted by a physician or specialist dealer so that actual pressure in the mask corresponds to therapy pressure.

4.3 Exhalation system

Vented variants

Vented variants (transparent retaining ring and elbow) have an integrated exhalation system. The retaining ring and mask body are shaped so that there is a gap between these parts. The exhaled air can escape through this gap.

Non-vented variants

Non-vented variants (blue retaining ring) do not have an exhalation system. Only use non-vented full-face masks in combination with therapy devices which have an active exhalation valve and have alarms and safety systems for any device failure. If you are using the non-vented variant with an external exhalation valve, please see the associated instructions for use.

4.4 Anti-asphyxia valve

⚠ WARNING

Risk of asphyxia if anti-asphyxia valve not working properly!

Residues may cause the valve to stick and lead to CO₂ being re-inhaled.

⇒ Check before every use that the openings of the anti-asphyxia valve are clear.

If the therapy device fails, the anti-asphyxia valve opens so that the patient can breathe ambient air.

5 Hygiene treatment

The mask is intended for a single patient and a maximum usage period of 7 days. A change of patient is not permitted.

Clean the mask with an isopropyl alcohol wipe (70 % by vol.) before using it for the first time and if it is lightly soiled.

6 Disposal

Dispose of the mask in accordance with hospital regulations.

7 Troubleshooting

Fault	Cause	Remedy
Pain from pressure on the face.	Mask too tight.	Loosen headgear slightly.
Draft in the eye.	Mask too loose.	Tighten headgear slightly.
	Mask does not fit.	Use mask of different size.
Therapy pressure is not reached.	Mask not correctly adjusted.	Re-adjust mask (see Figure 2).
	Patient circuit leaking.	Check connectors, check tubes properly located.

8 Technical data

Product class to directive 93/42/EEC	IIa
Dimensions (W x H x D)	
Size S	100 mm x 145 mm x 110 mm
Size M	100 mm x 155 mm x 110 mm
Size L	100 mm x 165 mm x 110 mm
Weight	
Size S	110 g
Size M	113 g
Size L	117 g
Therapy pressure range	4 hPa - 50 hPa
Tube connection: tapered connector to EN ISO 5356-1	Ø 22 mm (female)
Temperature range:	
Operation	+5 °C to +40 °C
Transport and storage	-20 °C to +70 °C
Flow resistance, vented	
at 50 l/min	0.04 hPa
at 100 l/min	0.19 hPa
Flow resistance, non-vented	
at 50 l/min	0.04 hPa
at 100 l/min	0.11 hPa
Flow resistance, non-vented + AAV	
at 50 l/min	0.05 hPa
at 100 l/min	0.19 hPa
Flow resistance, anti-asphyxia valve	
Inspiration at 50 l/min:	0.8 hPa
Exhalation at 50 l/min:	0.5 hPa

Switching pressure	
Anti-asphyxia valve	
• Open:	1.2 hPa
• Close:	1.4 hPa
Geometrical dead space:	
Size S	225 ml
Size M	255 ml
Size L	305 ml
Service life	5 years
Useful life	Up to 7 days

9 Materials

Only use the mask after consulting your physician if you are allergic to any of these substances:

Mask part	Material
Mask cushion, forehead cushion	SI (Silicone)
Mask body	HTC (Tritan)
Retaining ring	PP (polypropylene)
Elbow	PC (polycarbonate)
Headgear clip	POM (polyoxymethylene)
Quick-release cord	POM (polyoxymethylene), CO (cotton), PU (polyurethane)
Headgear	EL (elastane), PES (polyester), PU (polyurethane), UBL loop, CO (cotton), PA (polyamide)

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

10 Warranty

Löwenstein Medical Technology gives the customer a limited manufacturer warranty on new original Löwenstein Medical Technology products and any replacement part fitted by Löwenstein Medical Technology in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions are available on the website of the manufacturer. We will also send you the warranty conditions on request.

Please bear in mind that any claim to warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine replacement parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

11 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The unabbreviated text of the Declaration of Conformity can be found on the manufacturer's website.

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