
EN-US Instructions for use PL Instrukcja używania HU Használati utasítás CS Návod k použití SK Návod na použitie HR Upute za upotrebu RO Instrucțiuni de utilizare BG Инструкции за употреба SL Navodila za uporabo SR Uputstvo za upotrebu LV Lietošanas instrukcija ET Kasutusjuhend LT Naudojimo instrukcija RU Инструкция по пользованию MK Упатство за употреба TH คู่มือการใช้งาน MS Arahan penggunaan VI Hướng dẫn sử dụng ID Pedoman Penggunaan

WM 68056k 01/2026 EN-US, PL, HU, CS, SK, HR, RO, BG, SL, SR, LV, ET, LT, RU, MK, TH, MS, VI, ID

CE 0197

 **Manufacturer**  
Löwenstein Medical  
Technology GmbH + Co.KG  
Kronsaalweg 40  
22525 Hamburg, Germany  
T: +49 40 54702-0  
F: +49 40 54702-461  
[www.loewensteinmedical.com](http://www.loewensteinmedical.com)








**JOYCEone Full Face, JOYCEone Full  
Face NV**  
Full Face Mask

WM 68056k

**LÖWENSTEIN**  
medical

**LÖWENSTEIN**  
medical

Symbol	Description
	Use by date
	Keep out of sunlight
	Lot number
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

## 11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor original spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

## 12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website ([www.loewensteinmedical.com](http://www.loewensteinmedical.com)).

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

## 1 Operation

### Operating steps

The figures show the following steps for operating the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask

### For blind and partially-sighted users

An electronic version of the instructions for use is also available from the download area of the manufacturer's website.

## 2 Introduction

### 2.1 Intended use

The JOYCEone Full Face/JOYCEone Full Face NV mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device. Use non-vented full-face masks only in combination with therapy devices that have an active exhalation system.

### 2.2 Contraindications

The mask must not be used on patients weighing < 30 kg.

The mask must not be used in the following situations: Immediate intubation required, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: Pressure marks and acute injuries to the skin of the face; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea.

The mask is not suitable for use with a nebulizer or as an anesthetic mask. If you are not sure whether one of these situations applies to you, consult your healthcare professional. Observe the contraindications in the instructions for use of your device.

### 2.3 Side effects

The following side effects may occur when using the mask: nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the paranasal sinuses, irritated conjunctiva, skin redness, pressure marks on the face, facial edema, irritating noises when breathing.

Contact your healthcare professional if any of these side effects occurs.

### 2.4 Clinical benefit

Transfers the therapeutic efficacy of the therapy device to the patient.

## 3 Safety

**⚠ WARNING!** Unusually hazardous situation. If you ignore the following instructions, severe, irreversible or fatal injuries may result.

### Risk of injury due to damaged mask parts or those under strain!

- ⇒ Perform a visual inspection before every use and after every cleaning operation. Replace mask parts if necessary.
- ⇒ Note useful life (see section entitled "Technical data").

### Risk of injury from accumulation of CO<sub>2</sub>!

Expired carbon dioxide (CO<sub>2</sub>) may accumulate in the mask.

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.
- ⇒ Check before every use that the openings of the anti-asphyxia valve are clear.
- ⇒ Do not close off exhalation systems.

### Risk of injury due to the use of oxygen!

Oxygen can accumulate in clothing, bed linen, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ If required: Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.
- ⇒ Observe the safety instructions in the instructions for use of the device.

### Risk of injury due to insufficient supply to the patient!

- ⇒ Activate low pressure/leakage alarms on device.
- ⇒ Use the correct mask size and check that it fits tightly.
- ⇒ Monitor patients with restricted spontaneous breathing.

### Risk of injury from inadequate cleaning!

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.
- ⇒ In case of a change of patient in the clinical environment: Follow the document *Mask reprocessing instructions* (see section entitled "Reprocessing").

**⚠ CAUTION!** Hazardous situation. If you ignore these instructions, mild or moderate injuries may result.

### Risk of injury due to allergies!

- ⇒ Be aware of the materials used in the mask (see section entitled "Materials"). Only use a mask following consultation with a healthcare professional.
- ⇒ When selecting the detergent, consider potential allergies.

### Risk of injury due to incorrect adjustment of the headgear!

If the headgear is tightened too much, the mask may fit too tightly, press against the face, and cause skin redness, pressure marks, or facial edema.

- ⇒ Adjust the headgear correctly (see Figure 2) so that the mask fits comfortably and seals without excessive pressure.

## 4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	6	Anti-asphyxia valve
2	Forehead cushion	7	Rotating sleeve
3	Forehead support	8	Mask body
4	Retaining ring	9	Headgear clip
5	Elbow	10	Mask cushion

### Compatible devices

In some device combinations, the actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that the actual pressure in the mask corresponds to the therapy pressure. This adjustment should be made with the type of mask that will be used during therapy.

### Exhalation system

Masks with an integrated exhalation system have a gap through which exhaled air escapes.

Only use masks without an integrated exhalation system ("NV", blue retaining ring and elbow) with devices that have an active exhalation system as well as alarms and safety systems for the event of the device failing. If using external exhalation systems, follow the associated instructions for use.

### Anti-asphyxia valve (AAV)

If the device fails, the anti-asphyxia valve opens and the patient breathes ambient air.

### Quick-release cord (optional)

The quick-release cord allows the mask to be released quickly and easily in emergency situations (see SOS illustration).

## 5 Cleaning and reprocessing

### 5.1 Clean mask

1. Wash your hands before starting cleaning.
2. Dismantle mask (see Figure 4).

3. Clean the mask by hand (max. 30 °C, 1 ml mild detergent to 1 l water) in accordance with the following table:

Mask part	Frequency	Action
All mask parts	Daily	Soak for 15 minutes and wash for 3 minutes. Clean for 3 minutes with a soft cleaning brush.
Headgear	Weekly	Wash for 15 minutes.

All parts (exception: elbow with anti-asphyxia valve) can be washed in a dishwasher once a week (max. 70 °C, mild dishwashing detergent, max. 90 minutes program duration, top basket, separate rinse).

4. Rinse all parts with clean water.
5. Allow all parts to air-dry.
6. Carry out a visual inspection for cracks and deformations. Replace damaged parts. Discoloration is harmless.

7. Re-assemble mask (see Figure 5).

### 5.2 Reprocessing (clinical environment)

In the event of a change of patient, follow the instructions in the document entitled *Information on reprocessing*. The document can be found on the manufacturer's website. We will send you this document on request.

### 5.3 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

## 6 Troubleshooting

Fault	Cause	Action
Draft in the eye	Mask fits too loosely.	Adjust headgear so that it is tighter.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask is not correctly adjusted.	Re-adjust mask.
	Mask cushion is damaged.	Replace mask cushion.
	Circuit is damaged.	Check circuit and the correct fit of circuit.
Therapy pressure not reached.	Anti-asphyxia valve defective.	Replace anti-asphyxia valve.

## 7 Technical data

	Vented	NV
Dimensions in mm (H x W x D)	89 x 155 x 97	89 x 155 x 97

	Vented	NV
Weight	127 g	127 g
Dead space volume	240 ml	240 ml
Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (male)	Ø 22 mm (female)
Flow resistance at 50 l/min at 100 l/min	0.3 hPa 0.7 hPa	0.03 hPa 0.11 hPa
AAV flow resistance Insp. at 50 l/min Exp. at 50 l/min Tolerance: ± 0.2 hPa	0.6 hPa 0.8 hPa	- -
AAV switching pressure Open Close	≤ 1.0 hPa < 4.0 hPa	- -

Service life	5 years
Useful life	Up to 12 months <sup>1</sup>
Therapy pressure	4 hPa - 30 hPa
Quoted two-figure noise emission value according to ISO 4871: Sound pressure level Sound power level Uncertainty factor	19 dB(A) 27 dB(A) 3 dB(A)
Temperature: Operation Transport and storage	+5 °C to +40 °C -20 °C to +70 °C
Standards applied	EN ISO 17510: 2020
Product class to MDR (EU) 2017/745	IIa

<sup>1</sup> Useful life depends on cleaning and on the cleaning agents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. The use of the second mask cushion (optionally included) will not extend the useful life.


## 8 Materials

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

Headgear	CO (cotton), EL (elastane), PA (polyamide), P (polyester), PU (polyurethane)
Headgear clip	PA (polyamide), POM (polyoxymethylene)
Quick-release cord	PET (polyethylene terephthalate fiber), PA (polyamide)

Quick-release cord clip	PA (polyamide), POM (polyoxymethylene)
Forehead support	PA (polyamide)
Spring of forehead support	SI (silicone)
Forehead cushion	SI (silicone)
Retaining ring (hard component)	PP (polypropylene)
Retaining ring (soft component)	TPE (thermoplastic elastomer)
Mask body	PA (polyamide)
Mask cushion	SI (silicone)
Elbow	PA (polyamide)
Rotating sleeve	PA (polyamide)
Anti-asphyxia valve	SI (silicone)
Valve safety device	PP (polypropylene)

## 9 Characteristic pressure/flow curve

## 10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Unique device identifier (unique product code for medical devices)
	Order number
	Indicates the product is a medical device
	Follow the instructions for use
	Permitted temperature range for transport and storage

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