
DE Gebrauchsanweisung EN-US Instructions for use FR Mode d'emploi IT Istruzioni d'uso
ES-MX Manual de instrucciones PT Instruções de uso NL Gebruiksaanwijzing DA Brugsan-
visning NO Bruksanvisning SV Bruksanvisning FI Käyttöohje TR Kullanma talimatı EL
Οδηγίες χρήσης KO 사용 설명서 AR تعليمات الاستخدام

LMT 67690d 03/2026 DE, EN-US, FR, IT, ES-MX, PT, NL, DA, NO, SV,
FI, TR, EL, KO, AR

Printed on 100 % recycled paper

CE 0197

 **Manufacturer**
Löwenstein Medical
Technology GmbH + Co.KG
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.com







JULIA
Nasal Mask


LMT 67690d

LÖWENSTEIN
medical

LÖWENSTEIN
medical

9 Druck-Flow-Kennlinie

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Symbol	Description
	Order number
	Indicates the product is a medical device
	Consult instructions for use
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Lot number
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please be aware that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor original spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website (www.loewensteinmedical.com).

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

1 Operation

Before first use

Clean the mask before first use (see chapter Cleaning and reprocessing).

Steps for operation

The following steps for operation of the mask are shown in the figures:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask

For blind and partially-sighted users

An electronic version of the instructions for use is also available from the download area of the manufacturer's website.

2 Introduction

2.1 Intended use

The nasal mask is used as an accessory for devices for treating sleep-related breathing disorders or for non-life-support ventilation in the case of respiratory insufficiency. The nasal mask is suitable for reuse for patients in the home environment. In the clinical environment, the nasal mask is suitable for multiple patients after reprocessing.

2.2 Contraindications

The mask must not be used on patients weighing < 30 kg.

The mask must not be used in the following situations: Immediate intubation required, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: Pressure marks and acute injuries to the skin of the face; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea. The mask is not suitable for use with a nebulizer or as an anesthetic mask.

If you are not sure whether one of these situations applies to you, consult your healthcare professional. Observe the contraindications in the instructions for use of your device.

2.3 Side effects

The following side effects may occur when using the mask: nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the paranasal sinuses, irritated conjunctiva, skin redness, pressure marks on the face, facial edema, irritating noises when breathing.

Contact your healthcare professional if any of these side effects occurs.

2.4 Clinical benefit

Transfers the therapeutic efficacy of the therapy device to the patient.

3 Safety

⚠ WARNING! Unusually hazardous situation. If you ignore the following instructions, severe, irreversible or fatal injuries may result.

Risk of injury due to damaged mask parts or those under strain!

- ⇒ Perform a visual inspection before every use and after every cleaning operation. Replace mask parts if necessary.
- ⇒ Note useful life (see section entitled "Technical data").

Risk of injury from accumulation of CO₂!

Expired carbon dioxide (CO₂) may accumulate in the mask.

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Do not use the mask on patients who are unable to take it off themselves.
- ⇒ Do not close off exhalation systems.

Risk of injury due to the use of oxygen!

Oxygen can accumulate in clothing, bed linen, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ If required: Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.
- ⇒ Observe the safety instructions in the instructions for use of the device.

Risk of injury due to insufficient supply to the patient!

- ⇒ Activate low pressure/leakage alarms on device.
- ⇒ Use the correct mask size and check that it fits tightly.

Risk of injury from inadequate cleaning!

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.
- ⇒ In case of a change of patient in the clinical environment: Follow the relevant reprocessing instructions to your mask (see section entitled "Reprocessing").

Risk of injury if tubes are routed incorrectly!

- ⇒ Do not route tubes along the neck.
- ⇒ Do not crush tubes.

⚠ CAUTION! Hazardous situation. If you ignore these instructions, mild or moderate injuries may result.

Risk of injury due to allergies!

- ⇒ Be aware of the materials used in the mask (see section entitled "Materials"). Only use a mask following consultation with a healthcare professional.
- ⇒ When selecting the detergent, consider potential allergies.

Risk of injury due to incorrect adjustment of the headgear!

If the headgear is tightened too much, the mask may fit too tightly, press against the face, and cause skin redness, pressure marks, or facial edema.

- ⇒ Adjust the headgear correctly (see Figure 2) so that the mask fits comfortably and seals without excessive pressure.

4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	6	Adapter (elbow to mask tube)
2	Frame (left) and frame (right), each with frame cover	7	Mask tube
3	Mask cushion (two mask cushions included, depending on variant)	8	Adapter (rotating sleeve to breathing tube)
4	Mask body	9	Rotating sleeve
5	Elbow		

Compatible devices

In some device combinations, the actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that the actual pressure in the mask corresponds to the therapy pressure. This adjustment should be made with the type of mask that will be used during therapy.

Exhalation system

The mask has an integrated exhalation system. The exhaled air escapes through a gap.

5 Cleaning and reprocessing

5.1 Clean mask

1. Wash your hands before starting cleaning.
2. Dismantle mask (see Figure 4).

i Leave the frame covers on the frames and clean them together. If they are very dirty or damaged, remove the frame covers and replace them.

3. Clean the mask by hand at a maximum water temperature of 30°C. Follow the manufacturer's instructions regarding concentration and exposure time of the cleaning agent.

Mask part	Frequency	Action
Headgear, frame with frame cover	Weekly	Wash for 15 minutes.
Mask body, elbow, mask cushion, rotating sleeve	Daily	Soak for 15 minutes and wash. Clean for 3 minutes with a soft cleaning brush.
Mask tube	Daily	Soak and wash for 15 minutes
	Weekly	After daily cleaning, clean the inside of the mask tube for 3 minutes with a soft cleaning brush.

i Mask body with elbow, mask cushion, and rotating sleeve can be cleaned weekly in the dishwasher (max. 70 °C, mild dishwashing detergent, max. 90 minutes program duration, upper basket, separate rinse cycle).

4. Rinse all parts with clean water.
5. Allow all parts to air-dry.
6. Hang up the mask tube to dry it.
7. Carry out a visual inspection for cracks and deformations. Replace damaged parts. Discoloration is harmless.
8. If the mask cushion is damaged or heavily soiled: Replace mask cushion. If a second mask cushion is not enclosed, contact your specialist dealer.
9. Re-assemble mask (see Figure 5).

5.2 Reprocessing (clinical environment)

In the event of a change of patient, follow the relevant reprocessing instructions for your mask. The document can be found on the manufacturer's website. We will send you this document on request.

5.3 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

6 Troubleshooting

Fault	Cause	Action
Draft in the eye	Mask fits too loosely.	Adjust headgear so that it is tighter.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask is not correctly adjusted.	Re-adjust mask.
	Mask cushion is damaged.	Replace mask cushion.
	Circuit is damaged.	Check circuit and the correct fit of circuit.

7 Technical data

Dimensions in mm (W x H x D)	
Size S	53 x 61 x 75
Size M	58 x 60 x 76
Size L	62 x 64 x 77
Weight	
Size S	77 g
Size M	78 g
Size L	80 g
Dead space	
Size S	68 ml
Size M	78 ml
Size L	92 ml
Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (male)
Flow resistance at 50 l/min at 100 l/min	0.6 hPa 2.3 hPa
Service life	5 years
Useful life	Up to 12 months ¹
Therapy pressure	4 hPa - 30 hPa
Quoted two-figure noise emission value according to ISO 4871:	
Sound pressure level	18 dB(A)
Sound power level	26 dB(A)
Uncertainty factor	3 dB(A)
Temperature: Operation	+5 °C to +40 °C
Transport and storage	-20 °C to +70 °C
Standards applied	EN ISO 17510: 2020
Product class to MDR (EU) 2017/745	IIa


¹ Useful life depends on cleaning and on the cleaning agents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. The use of the second mask cushion (optionally included) will not extend the useful life.

8 Materials

All parts of the mask are free from latex, PVC (polyvinyl chloride), DEHP (diethylhexyl phthalate) and BPA (bisphenol-A).

Headgear	PA (polyamide), PU (polyurethane), P (polyester)
Mask body	PA (polyamide)
Mask cushion	SI (silicone)
Mask tube	TPO (thermoplastic polyolefin elastomer)
Adapter (elbow to mask tube)	MABS (methyl methacrylate acrylonitrile butadiene styrene)
Adapter (rotating sleeve to breathing tube)	MABS (methyl methacrylate acrylonitrile butadiene styrene)
Elbow	PA (polyamide)
Rotating sleeve	PA (polyamide)
Frame (left/right)	PA (polyamide)
Frame cover	P (polyester)

9 Characteristic pressure/flow curve

10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Unique device identifier (unique product code for medical devices)

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