
EN-US Instructions for Use PT-BR Instruções de uso HR Uputa za uporabu EL Οδηγίες χρήσης PL Instrukcja używania CS Návod k použití SK Návod na použitie HU Használati utasítás RO Instrucțiuni de utilizare SR Uputstvo za upotrebu SL Navodila za uporabo MK Παττβο за употреба DA Brugsanvisning NO Bruksanvisning SV Bruksanvisning FI Käyttöohje KO 사용 설명서 MS Arahan penggunaan

LMT 65881e 12/2023 EN-US, PT-BR, HR, EL, PL, CS, SK, HU, RO, SR, SL, MK, DA, NO, SV, FI, KO, MS

CE 0197

 **Manufacturer**  
Löwenstein Medical  
Technology GmbH + Co.KG  
Kronsaalweg 40  
22525 Hamburg, Germany  
T: +49 40 54702-0  
F: +49 40 54702-461  
[www.loewensteinmedical.com](http://www.loewensteinmedical.com)









**LENA, LENA NV**  
Full Face Mask

LMT 65881e

**LÖWENSTEIN**  
medical

**LÖWENSTEIN**  
medical

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Follow the instructions for use
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Lot number

## 11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor genuine spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

## 12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

## 1 Operation

The figures show the following steps for operating the mask:

- 1** Putting on the mask
- 2** Adjusting the mask
- 3** Removing the mask
- 4** Dismantling the mask
- 5** Assembling the mask
- 6** Supplying oxygen

**i** **For blind and partially-sighted users**  
An electronic version of the instructions for use is also available on the manufacturer's website.

## 2 Introduction

### 2.1 Intended use

LENA is designed for use as an accessory by patients weighing over 30 kg who have been prescribed non-invasive positive airway pressure therapy (PAP therapy) - CPAP/APAP, BiLevel or NIV, for example. The mask is suitable for reuse for individual patients in the home environment and for reuse for individual and multiple patients in the hospital/care home environment.

### 2.2 Contraindications

The mask must not be used on patients weighing < 30 kg.

The mask must not be used in the following situations: Immediate intubation required, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: Pressure points and acute injuries to the skin of the face; skin allergies involving the face; deformities of the face or nasopharynx; acute pain affecting the face; cough reflex restricted or absent; claustrophobia; acute nausea.

If you are not sure whether one of these situations applies to you, consult your healthcare professional. Observe the contraindications in the instructions for use of your device.

### 2.3 Side effects

The following side effects may occur with use of the mask: Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

If these side effects occur, contact your healthcare professional.

### 2.4 Clinical benefit

Transfers the therapeutic efficacy of the ventilator to the patient

## 3 Safety

**Risk of injury due to damaged mask parts or those under strain!**

- ⇒ Perform a visual inspection before every use and after every cleaning operation.
- ⇒ Note useful life (see section entitled "Technical Data").
- ⇒ Replace mask parts if necessary.

**Risk of injury due to the use of oxygen!**

Oxygen can become deposited in clothing, bedding, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Set up oxygen sources at a distance > 1 m from the device.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.

**Risk of injury due to patient receiving inadequate supply!**

- ⇒ Activate low pressure/leakage alarms on the ventilator.
- ⇒ Use the appropriate mask size and check that it is securely in position.
- ⇒ Monitor patients with restricted spontaneous breathing.

**Risk of injury from re-inhalation of CO<sub>2</sub>!**

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients unable to remove the mask themselves must be monitored by a nurse.
- ⇒ Check before every use that the openings of the anti-asphyxia valve are clear.
- ⇒ Do not close off exhalation systems.

**Risk of injury due to escape of anesthetic gas or atomization of drugs!**

- ⇒ Do not use the mask during anesthesia.
- ⇒ Do not use the mask to atomize drugs.

**Risk of injury from inadequate cleaning!**

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.
- ⇒ When selecting a detergent, consider potential allergies.
- ⇒ On change of patient in a hospital environment: Comply with the document entitled "Information on reprocessing" (see section entitled "Reprocessing").
- ⇒ For patients with a compromised immune system or particular background of illness, disinfect mask parts daily following consultation with the healthcare professional.

## 4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	7	Elbow
2	Forehead cushion	8	Anti-asphyxia valve
3	Forehead support	9	Rotating sleeve
4	Adjusting element	10	Mask body
5	O <sub>2</sub> inlet	11	Headgear clip
6	Retaining ring	12	Mask cushion (two masks may be included, depending on variant)

### Compatible devices

In some device combinations, actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that actual pressure in the mask corresponds to therapy pressure. This setting should be made using the type of mask used during therapy.

### Exhalation system

Masks with an integrated exhalation system have a gap through which exhaled air escapes.

Only use masks without an integrated exhalation system ("NV", blue retaining ring and elbow) with devices that have an active exhalation system as well as alarms and safety systems for the event of the device failing. If using external exhalation systems, follow the associated instructions for use.

### Anti-asphyxia valve (AAV)

If the device fails, the anti-asphyxia valve opens and the patient breathes ambient air.

### Quick-release cord (optional)

The quick-release cord allows the mask to be released quickly and easily in emergency situations (see SOS illustration).

## 5 Cleaning and reprocessing

### 5.1 Clean mask

1. Wash your hands before starting cleaning.
2. Dismantle mask (see Figure 4).
3. Wash mask by hand (max. 30 °C, 1 ml mild detergent in 1 l water) in accordance with the table below:

Mask part	Frequency	Action
All mask parts	Daily	Soak and wash for 15 minutes and clean for 3 minutes using a soft cleaning brush.
Headgear	Weekly	Wash for 15 minutes.

**i** All parts (exception: elbow and anti-asphyxia valve) can be washed in a dishwasher once a week (max. 70 °C, mild detergent, max. program length 90 minutes, top basket, separate rinse).

4. Rinse all parts with clean water.
5. Allow all parts to air-dry.
6. Perform a visual inspection for cracks and deformation. Replace damaged parts. Discoloration is not a problem.
7. If the mask cushion is damaged or heavily soiled: Replace mask cushion. If a second mask cushion is not enclosed, contact your specialist dealer.
8. Assemble mask (see Figure 5).

### 5.2 Reprocessing (clinical environment)

In the event of a change of patient, follow the instructions in the document entitled *Information on reprocessing*. The document can be found on the manufacturer's website. We will send you the document on request.

### 5.3 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

## 6 Troubleshooting

Fault	Cause	Action
Pain due to pressure on the face	Mask too tight.	Loosen headgear.
Draft in the eye	Mask too loose.	Tighten headgear.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure not reached.	Mask not correctly adjusted.	Re-adjust mask.
	Mask cushion damaged.	Replace mask cushion.
	Patient circuit damaged.	Check circuit and correct fit of circuit.
Therapy pressure not reached.	Anti-asphyxia valve defective.	Replace anti-asphyxia valve.

## 7 Technical data

	Vented	NV
Dimensions in mm (H x W x D) Size S Size M Size L	155 x 100 x 95 165 x 100 x 95 175 x 100 x 100	155 x 100 x 105 165 x 100 x 105 175 x 100 x 110
Weight Size S Size M Size L	137 g 141 g 150 g	135 g 139 g 148 g
Dead space Size S Size M Size L	246 ml 288 ml 326 ml	252 ml 270 ml 321 ml
Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (male)	Ø 22 mm (female)
Flow resistance at 50 l/min at 100 l/min	0.32 hPa 0.67 hPa	0.04 hPa 0.14 hPa
AAV flow resistance Insp. at 50 l/min Exp. at 50 l/min Tolerance: ± 1 hPa	0.6 hPa 0.8 hPa	- -
AAV switching pressure Open Close	0.5 hPa 2.2 hPa	- -


Service life	5 years
Useful life	Up to 12 months <sup>1</sup>
Therapy pressure	4 hPa - 35 hPa
Quoted two-figure noise emission value according to ISO 4871: Sound pressure level Sound power level Uncertainty factor	12 dB(A) 20 dB(A) 3 dB(A)
Temperature: Operation Transport and storage	+5 °C to +40 °C -20 °C to +70 °C
Standards applied	EN ISO 17510: 2020
Product class to MDR (EU) 2017/745	IIa
<sup>1</sup> Useful life depends on cleaning and on the detergents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. Using the second mask cushion (included as an option) does not extend useful life.	

## 8 Materials

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate). In the event of allergies to the materials listed, only use the mask following agreement with the healthcare professional.

Headgear	PA (polyamide), P (polyester), PU (polyurethane)
Headgear clip	POM (polyoxymethylene)
Quick-release cord	PET (polyethylene terephthalate fiber), PA (polyamide)
Quick-release cord clip	POM (polyoxymethylene)
Adjusting element	POM (polyoxymethylene)
O <sub>2</sub> inlet	SI (silicone)
Retaining ring	POM (polyoxymethylene)
Mask body	PA (polyamide)
Mask cushion	SI (silicone)
Elbow	PA (polyamide)
Rotating sleeve	PP (polypropylene)
Anti-asphyxia valve	SI (silicone)
Valve safety device	PP (polypropylene)

## 9 Characteristic pressure/flow curve

## 10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
<b>UDI</b>	Unique device identifier (uniform product code for medical devices)
<b>REF</b>	Order number
<b>MD</b>	Indicates the product is a medical device

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