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
10.5Warranty 53
10.6Declaration of Conformity 54

1 Intended use

The COM cable for monitor is for transmitting data from LM150TD ventilators to Philips IntelliVue monitors.

 LM150TD ventilators include LUISA, TIVAN LS and LifeOne.

2 Safety

2.1 Safety information

2.1.1 Handling the product and its components

If the product is damaged or its function is restricted, individuals may be injured.

- ⇒ The use of third-party articles may lead to incompatibility with the device. In such cases, please be aware that any claim under warranty and liability will be void if original spare parts are not used.
- ⇒ Connect only the licensed products and modules in accordance with these instructions for use. The products must meet the product standard applicable to them. Non-medical equipment should be positioned out of the patient's vicinity.
- ⇒ Only use the product if the exterior is undamaged.
- ⇒ Perform a function check at regular intervals (see ["4.2 Function check"](#), page 37).
- ⇒ Only use, store and transport the product under the specified ambient conditions (see ["10.1 Technical specifications"](#), page 52).
- ⇒ Keep small parts which may be inhaled or swallowed away from young children in particular.
- ⇒ Do not use the product in an MRI environment or in a hyperbaric chamber.

2.1.2 Electromagnetic compatibility

- ⇒ The product is subject to special precautions with regard to EMC (electromagnetic compatibility). If these precautions are not followed, the product may malfunction and individuals may be injured.
- ⇒ Portable high-frequency communication equipment (e.g. radios and cell phones), including their accessories such as antenna cables and external antennas, for example, must be used at a distance of at least 30 cm from the product.
- ⇒ Do not use the product in the vicinity of active high-frequency surgical equipment.
- ⇒ Do not operate the product if the housing, cables or other equipment for electromagnetic shielding are damaged.

- ⇒ Do not operate the product in the immediate vicinity of other products or in a stacked arrangement, otherwise there may be malfunctions. If it is necessary to operate the product in the immediate vicinity of other products or in a stacked arrangement, keep all the products under observation to ensure that they are all operating properly.
- ⇒ For more information on EMC, follow the instructions for use for the ventilator.

2.2 General information


- The COM cable for monitor is intended for clinical use.
- The instructions for use for the device, the components, and the accessories should also be followed.
- Have measures such as repairs, servicing, and maintenance work, as well as modifications to the product carried out exclusively by the manufacturer or by specialists expressly so authorized by the manufacturer.
- Follow the section on hygiene treatment (see "6 Hygiene treatment", page 50) to avoid infection or bacterial contamination.
- In the EU: As a user and/or patient, you must report any serious incident occurring in conjunction with the product to the manufacturer and to the responsible authority.

2.3 Warnings in these instructions for use

 WARNING	Indicates an unusually significant hazardous situation. If you ignore this instruction, severe irreversible or fatal injuries may result.
 CAUTION	Indicates a hazardous situation. If you ignore this instruction, mild or moderate injuries may result.
NOTICE	Indicates a harmful situation. If you ignore this instruction, material damage may result.
	Indicates useful information within procedures.

3 Product description

3.1 Overview

1. COM cable for monitor LMT 31578
2. LM150TD ventilator
3. Philips IntelliBridge EC5 module, including connecting cable [article number Philips 865114/101/L01 (1.5 m), 865114/101/L02 (3 m), 865114/101/L03 (10 m)]
4. Philips IntelliVue monitor
5. Philips IntelliBridge EC10 module (article number Philips 865115, option A01/101)

3.2 Compatibility

The LM150TD ventilator is compatible with:

Philips Patient Monitoring Systems IntelliVue MP40-90 and MX400-850 with Philips IntelliBridge EC10 Interface Module or IntelliBridge EC10 integral Interface Board

4 Commissioning

4.1 Operation

⚠ CAUTION

Risk of injury if cables are routed incorrectly!

Incorrectly routed cables may injure the patient.

⇒ Do not route cables along the neck and keep them away from the patient.

NOTICE


Risk of material damage!

Incorrect use may damage the product.


⇒ Only use the COM cable for monitor to connect the device to the monitor.

Requirement


- The device and the patient monitor are switched on.
1. Insert the IntelliBridge EC10 module in the patient monitor (see associated instructions for use).

2. Connect the COM cable for monitor to the device.
3. Connect the IntelliBridge EC5 connecting cable to the gray IntelliBridge EC5 adapter.

4. Connect the IntelliBridge EC5 connecting cable to the IntelliBridge EC10 module.

5. Connect the IntelliBridge EC5 connecting cable to the COM cable for monitor.
The connection from the device to the patient monitor is complete. It takes about 45 seconds to set up the connection.

4.2 Function check

Carry out a function check before using the device for the first time, after every hygiene treatment, and after every repair, but at least every 6 months.

Performing a visual inspection

1. Check product for external damage.

Checking data transfer

1. Connect the device to the patient monitor (see "4.1 Operation", page 35).
2. Start therapy on the device.
3. Requirement: Values or signals from the device will be displayed on the patient monitor.

Checking the disconnection alarm

1. Uncouple the connection between the device and the patient monitor.
2. Requirement: The device display will show the alarm **Disconnected from patient monitor** (Error number 788).

Checking alarms

1. Connect the device to the patient monitor.
2. Generate any alarm on the device (see the instructions for use for the device).

Example alarm *581 Power supply via internal battery*:

Disconnect the power cord and external batteries from the device.

3. Requirement: Alarm is displayed on the patient monitor.

If one of the items is not OK: Do not use the product and contact your specialist dealer.

5 Data mapping

These comparisons show how data elements are visualized in the two different data models, i.e. how the data are displayed in the LM150TD device and on the patient monitor.

The specified name of a value or alarm in column *Intellivue name* is the English name that is displayed on the Philips patient monitor. If the Philips patient monitor is set to a different language, other Intellivue names may be displayed.

5.1 Therapy parameter

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Patient's current tidal volume 0 ml - 3,000 ml	All	Curve	0002-f0df	AWV	Airway volume wave	Curve
Patient's current respiratory flow rate -120 l/min to 250 l/min	All HFT: Device-Out-Flow is displayed	Curve	0002-50d4	AWF	Airway flow wave	Curve
Current actual pressure at patient 0 hPa to 100 hPa	All	Curve	0002-50f0	AWP	Airway pressure wave	Curve
Current pulse (if available) 0 bpm to 255 bpm	all - if device and sensor connected	Measurement	0002-480a	PULSE	Pulse Rate	Numerical
Current SpO ₂ (if available) 0 % to 100 %	all - if device and sensor connected	Measurement	0002-4bb8	SpO_2	Percent Oxyhemoglobin in Saturation - Arterial	Numerical
Current therapy duration 0 hr to 999 h	All	Measurement	0002-fc3c	VentTi	Ventilation Time	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Current leakage 0 l/min to 200 l/min	CPAP, S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE (Leakage and double circuit. On single patient valve ventilation: = device peak flow.)	Measurement	0002-5370	Leak	Leakage	Numerical
Current CPAP 4 hPa to 100 hPa	CPAP	Measurement	0002-50f4	CPAP	CPAP	Numerical
Current Te 0 sec to 60 sec	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Measurement	0002-f8a1	ExpTi	Expiratory Time	Numerical
Current I:E ratio 0 to 10	CPAP (only valve circuit), S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Measurement	0002-5118	I:E 1:	Inspired:Expired Ratio	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Current Ti 0 sec to 10 sec	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Measurement	0002-f8a3	InsTi	Inspiratory Time	Numerical
Patient's minute volume 0 l/min to 100 l/min	S, ST, autoST, T, PSV, aPCV, PCV, P-SIMV, V-SIMV, TTV_VAPS_AE (only leakage circuit)	Measurement	0002-5148	MV	Minute Volume	Numerical
Current respiratory frequency 0 rpm to 100 rpm	CPAP (only valve circuit), S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Measurement	0002-5022	RRaw	Airway Respiration Rate	Numerical
Current IPAP max 0 hPa to 100 hPa	all CPAP (only valve circuit)	Measurement	0002-5109	PIP	Peak Inspiratory Pressure	Numerical
Current EPAP/ PEEP 0 hPa to 100 hPa	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Measurement	0002-51a8	PEEP	Positive Endexpiratory Pressure	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Setting parameter for Ti/ Ti min 0.2 s to 4 s	Ti for T, aPCV, PCV, aVCV, VCV, MPVp. Ti min for S, ST, autoST, PSV, P-SIMV, V-SIMV, TTV-VAPS-AE	Setting	0402-f941	sInsTi	Setting: Inspiratory Time	Numerical
Setting parameter for target volume 30 ml to 3,000 ml	S, ST, autoST, T, PSV, aPCV, PCV, TTV_VAPS_AE (if VTC active), aVCV, VCV, V-SIMV, MPVv	Setting	0402-f98f	sTVsup	Setting: Target Support Volume	Numerical
Setting parameter for trigger on inspiration 1 to 10, (auto: cannot be displayed on the monitor)	CPAP (only valve circuit), S, ST, autoST, PSV, aPCV, aVCV, P-SIMV, V-SIMV, TTV_VA PS_AE (if Trigger = manual), MPVp, MPVv	Setting	0402-f99a	sTrig	Setting: Trigger Sensitivity	Numerical
Setting parameter for trigger on expiration 5 % to 95 %	CPAP (only valve circuit), S, ST, autoST, PSV, P-SIMV, V-SIMV, TTV_VAPS_AE	Setting	0402-fbb2	sFISec	Setting: Cycle Sensitivity Sets flow sensitivity for switch from inspiration to expiration	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Setting parameter for Ti max 0.2 s to 4 s	S, ST, autoST, PSV, P-SIMV, V-SIMV, TTV_VAPS_AE	Setting	0402-f9e0	sTimax	Setting: Maximum Inspiratory Time	Numerical
Setting parameter for F (background frequency) 0 rpm to 80 rpm	ST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, on autoST, TTV_VAPS_AE - if autoF = Off	Setting	0402-f9cc	sRRbak	Setting: Backup Respiration Rate of ventilator	Numerical
Setting parameter for pressure rise 1 to 4	S, ST, autoST, T, PSV, aPCV, PCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Setting	0402-f912	sPincR	Setting: Pressure Increase Rate	Numerical
Setting parameter for EPAP/PEEP 0 hPa to 25 hPa	S, ST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV (EPAP displayed) autoST, TTV_VAPS-AE (EPAPmin displayed)	Setting	0402-51a8	sPEEP	Setting: PEEP/CPAP	Numerical
Setting parameter IPAP 4 hPa to 60 hPa	S, ST, T, PSV, aPCV, PCV, P-SIMV, MPVp, MPVv	Setting	0402-fbb7	sIPAP	Setting: Inspiratory positive airway pressure	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Setting parameter for flow 5 l/min to 60 l/min	HFT	Setting	0402-f91b	sFlow	Setting: Flow	Numerical
Setting parameter for CPAP 4 hPa to 20 hPa	CPAP	Setting	0402-50f4	sCPAP	Setting: Continuous Positive Airway Pressure Value	Numerical
Active mode	all	Setting	0402-f921	svMode	Enumeration Setting: Ventilation Mode	Numerical
Device type	LUISA	Setting	0402-f997	sVent	Enumeration Setting: Ventilator Type	Numerical
Setting parameter for pressure support 4 hPa - 60 hPa	P-SIMV, V-SIMV	Setting	0402-f8bc	sPSV	Setting: Pressure Support Ventilation	Numerical
Alarm paused	all	Setting	0402-f8d8	sSilnc	Setting: Alarm Silence Status	Numerical
Inspiratory oxygen fraction 0 % to 100 %	all (if sensor present)	Measurement	0002-7498	FIO_2	Fractional Inspired Oxygen	Numerical
Tidal volume on exhalation (only double circuit) 0 ml to 3000 ml	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV	Measurement	0002-f0e1	TVex	Expired Tidal Volume	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Tidal volume on inspiration 0 ml to 3000 ml (on leakage circuit: general TV)	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Measurement	0002-f0e0	TVin	Inspired Tidal Volume	Numerical
Minute volume on exhalation (only double circuit) 0 l/min to 99 l/min	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV	Measurement	0002-514c	MVexp	Expired Minute Volume	Numerical
Minute volume on inspiration 0 l/min to 99 l/min	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV (only for valve circuit), MPVp, MPVv	Measurement	0002-5150	MVin	Inspired Minute Volume	Numerical
Apnea duration 0 sec to 1800 sec	MPVp, MPVv	Measurement	0002-5130	ApneaD	Apnea Duration	Numerical
Triggered ventilation ratio (0 % to 100 %)	ST, autoST, PSV, aPCV, aVCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Measurement	0002-fb3e	%Bsp:t	Ratio of spontaneous/patient triggered breathing	Numerical
IPAP (0 hPa to 100 hPa)	S, ST, autoST, T, PSV, aPCV, PCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Measurement	0002-5368	Peinsp	End Inspiratory Pressure	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Alarm threshold for pressure high 6 hPa to 65 hPa	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Alarm limit	-	PIP_LIM	Limit Peak Inspiratory Pressure	Numerical
Alarm thresholds for leakage 5 l/min to 200 l/min	CPAP, S, ST, autoST, T, PSV, aPCV, PCV, TTV_VAPS_AE	Alarm limit	-	Leak_LIM	Limit Leakage	Numerical
Alarm thresholds for respiratory frequency 1 rpm to 80 rpm	S, ST, autoST, PSV, aPCV, aVCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Alarm limit	-	RRaw_LIM	Limit Airway Respiration Rate	Numerical
Alarm thresholds for SpO ₂ 85 % to 100 %	All - if device and sensor connected	Alarm limit	-	SpO2_LIM	Limit Percent Oxyhemoglobin Saturation - Arterial	Numerical
Alarm thresholds for pulse 30 bpm to 250 bpm	All - if device and sensor connected	Alarm limit	-	PULSE_LIM	Limit Pulse Rate	Numerical
Alarm threshold for tidal volume on exhalation 0 ml to 3000 ml	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV (only double circuit)	Alarm limit	-	TVex_LIM	Limit Expired Tidal Volume	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Alarm threshold for tidal volume on inspiration 0 ml to 3000 ml	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, MPVp, MPVv, TTV_VAPS_AE	Alarm limit	-	TVin_LIM	Limit Inspired Tidal Volume	Numerical
Alarm threshold for high PEEP 4 hPa to 30 hPa	S, ST, autoST, T, PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV, TTV_VAPS_AE	Alarm limit	-	PEEP_LIM	Limit Positive End-Expiratory Pressure	Numerical
Alarm threshold for minute volume 0.1 l/min to 40 l/min	S, ST, autoST, T, PSV, aPCV, PCV, P-SIMV, V-SIMV, TTV_VAPS-AE, (only leakage circuit)	Alarm limit	-	MV_LIM	Limit Minute Volume	Numerical
Alarm threshold for minute volume on exhalation 0.1 l/min to 40 l/min	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV (only double circuit)	Alarm limit	-	MVexp_LIM	Limit Expired Minute Volume	Numerical
Alarm threshold for minute volume on inspiration 0.1 l/min to 40 l/min	PSV, aPCV, PCV, aVCV, VCV, P-SIMV, V-SIMV (only valve circuit)	Alarm limit	-	MVin_LIM	Limit Inspired Minute Volume	Numerical
Alarm threshold for FIO ₂ 20 % to 100 %	all (if sensor connected)	Alarm limit	-	FIO_2_LIM	Limit Fractional Inspired Oxygen	Numerical

LM150TD			Intellivue			
Device description	Therapy mode	Data type	MDIL text ID	Intellivue name	Intellivue description	Data type
Alarm threshold for apnea duration 0 sec to 1800 sec	CPAP (only leakage circuit), S, MPVv, MPVp	Alarm limit	-	ApneaD_LIM	Limit Apnea Duration	Numerical

5.2 Alarm messages

Alarms are displayed on the patient monitor.

Fault no.	ID	Group	Priority level	Intellivue enum.	Intellivue classification
450	Tidal volume low	Physiological alarm	High	TV LOW	Red
451	Tidal volume high	Physiological alarm	High	TV HIGH	Red
452	Frequency low	Physiological alarm	Medium	FREQ LOW	Yellow
453	Frequency high	Physiological alarm	Medium	FREQ HIGH	Yellow
454	Minute volume low	Physiological alarm	High	MV LOW	Red
455	Minute volume high	Physiological alarm	Medium	MV HIGH	Yellow
456	Pressure high	Physiological alarm	High	PRESSURE HIGH	Red
457	Pressure low	Physiological alarm	Medium	PRESSURE LOW	Yellow
458	Apnea	Physiological alarm	High	APNEA	Red
459	Leakage high	Physiological alarm	Medium	LEAK HIGH	Yellow
460	Disconnection device outlet port	Technical alarm	High	DISCONNECT	Red
461	Disconnection airway pressure	Technical alarm	High	DISCONNECT	Red

Fault no.	ID	Group	Priority level	Intellivue enum.	Intellivue classification
463	Disconnection exhalation module	Technical alarm	High	DISCONNECT	Red
464	Disconnection patient	Technical alarm	High	DISCONNECT	Red
465	Disconnection patient HFT mode	Technical alarm	Low	DISCONNECT	Red
469	PEEP high	Physiological alarm	High	PEEP HIGH	Red
470	Exhalation tidal volume low	Physiological alarm	High	TV LOW	Red
471	Exhalation tidal volume high	Physiological alarm	High	TV HIGH	Red
472	Exhalation minute volume low	Physiological alarm	High	MV LOW	Red
473	Exhalation minute volume high	Physiological alarm	Medium	MV HIGH	Yellow
474	Inspiration tidal volume low	Physiological alarm	High	TV LOW	Red
475	Inspiration tidal volume high	Physiological alarm	High	TV HIGH	Red
476	Inspiration minute volume low	Physiological alarm	High	MV LOW	Red
477	Inspiration minute volume high	Physiological alarm	Medium	MV HIGH	Yellow
490	SpO ₂ low	Physiological alarm	High	SPO2 LOW	Red
491	SpO ₂ high	Physiological alarm	Medium	SPO2 HIGH	Yellow
492	Pulse low	Physiological alarm	High	PULSE LOW	Red
493	Pulse high	Physiological alarm	Medium	PULSE HIGH	Yellow
494	FiO ₂ low	Physiological alarm	High	FIO2 LOW	Red
495	FiO ₂ high	Physiological alarm	High	FIO2 HIGH	Red

Fault no.	ID	Group	Priority level	Intellivue enum.	Intellivue classification
794	Therapy ended	Technical alarm	High	VENTILATION OFF	Red
798	Disconnection valve control pressure	Technical alarm	High	DISCONNECT	Red
Others	Other alarms	Technical alarm	High	GEN TECH ALARM	Red

All technical alarms on the device that do not have a direct equivalent on the patient monitor are displayed on the patient monitor as **GEN TECH ALARM**.

⇒ For specific information on the respective alarm, refer to the alarm message on the device.

6 Hygiene treatment

CAUTION

Risk of injury from electric shock!

Ingress of liquids may lead to a short-circuit, injure the user and damage the device.

- ⇒ Disconnect the product from the power supply before the hygiene treatment.
- ⇒ Do not immerse the product and its components in liquids.
- ⇒ Do not pour liquid over the product and its components.

6.1 General information

- Wear appropriate safety gear (e.g. safety gloves) for the disinfecting process.
- Follow the instructions for use for the disinfectant used. Solutions containing alcohol (25 g ethanol (94 %-strength), 35 g propan-1-ol per 100 g) are suitable. Recommended: Mikrozyd AF liquid or perform advanced Alcohol EP.

6.2 Intervals

INTERVAL	ACTION
Weekly	Clean and disinfect the product

6.3 Cleaning the COM cable for monitor

1. Disconnect the cable connections.
2. Disconnect the COM cable for monitor from the device.
3. Dust off the COM cable for monitor with a soft cloth.
4. If necessary: Wipe over the product with a damp cloth. Use water or a mild detergent solution. Alternatively, use the disinfectant.
5. Perform function check ([see "4.2 Function check", page 37](#)).

6.4 Disinfecting the COM cable for monitor

Disinfect the COM cable for monitor by wiping it with mikrozyd AF wipes, for example.

7 Faults

If you are unable to remedy faults with the aid of the table, or in the event of unexpected operation or an incident, contact your authorized specialist dealer to have the product replaced. To avoid exacerbating the damage, do not continue operating the product.

Fault	Cause	Remedy
No connection to device.	Device firmware version is out of date and incompatible.	Check firmware version and have it updated by your specialist dealer if necessary.
	Cables plugged together incorrectly.	Maintain correct sequence of cables and connect them correctly (see 4.1, p. 35).

8 Servicing

The COM cable for monitor requires no servicing. The product is designed for a service life of 10 years. If used in accordance with the intended use, the device requires no servicing during this period.

9 Transport, storage and disposal

9.1 Transport and storage

Store and transport the COM cable for monitor under the specified ambient conditions (see "10.1 Technical specifications", page 52).

9.2 Disposal

Do not dispose of the product in domestic waste. To dispose of it properly, contact a licensed, certified electronic waste disposal merchant. This address is available from your Environment Officer or from your local authority. The device packaging (cardboard and inserts) can be disposed of in paper recycling facilities.

10 Annex

10.1 Technical specifications

SPECIFICATION	PRODUCT
Product class to MDR (EU) 2017/745	I
Length	180 cm
Weight	91 g
Temperature range Operation Transport and storage	+5 °C to +40 °C -25 °C to +70 °C
Permitted humidity for operation, transport and storage	15 % to 93 % relative humidity, no condensation At temperatures above 35° C, a maximum relative humidity corresponding to a partial water vapor pressure of 50 hPa.
Air pressure range	700 hPa to 1,060 hPa, corresponds to an altitude of 3,000 m above MSL
Electrical rating	3.3 V DC
Power capacity	< 0.5 W
Electromagnetic compatibility (EMC)	Test parameters and limit values can be obtained from the manufacturer on request.

*Data apply only in combination with the licensed device.
The right to make design modifications is reserved.

10.2 Markings and symbols

The following symbols may be applied to the device, the device ID plate, accessories, or packaging.

SYMBOL	DESCRIPTION
	Order number
	Manufacturer and date of manufacture
IP22	Degree of protection against contact with a finger. Protection against vertically falling water drops when enclosure tilted up to 15°.

SYMBOL	DESCRIPTION
	Consult instructions for use
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)
	Do not dispose of the device in domestic waste
	Permitted humidity range for transport and storage
	Permitted temperature range for transport and storage
	Serial number
	Unique device identifier (uniform product code for medical devices)

10.3 Scope of delivery

PART	ARTICLE NUMBER
COM cable for monitor	LMT 31578
Instructions for use	LMT 68580

10.4 Accessories and spare parts

PART	ARTICLE NUMBER
COM cable for monitor	LMT 31578

10.5 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

In the event of a claim under warranty, contact your specialist dealer.

PRODUCT	WARRANTY PERIODS
Devices including accessories (Exception: Masks)	2 years
Masks including accessories, rechargeable batteries, batteries (unless quoted differently in the technical documentation), sensors, patient circuits	6 months
Disposable products	None

10.6 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

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