DE

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Legende

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


1 Description of equipment

Legend

1. Sleeve with 22.7 mm connection
2. Porous filter
3. Connecting piece with 22 mm connection for EN ISO 5356-1
4. Plug

Intended purpose

The muffling system is used in combination with a nasal mask as an accessory in CPAP- and Bi-level therapy.

The muffling system enables you to exhale quietly during treatment. The muffling system prevents the exhaled air from being inhaled again.

Description of function

The muffling system connects the air hose of the treatment device to the nasal mask.

Inhaled air flows through the mask/tube system. Exhaled air escapes via porous filter **2**.

2 Safety Requirements

Special markings on the device

Storage temperature

Safety Instructions

For your own protection and that of your patient and according to the requirements of the European Union Guidelines 93/42/EEC, please observe the following:

- Please read these instructions for use carefully. They are a part of the device and have to be available at all times.
- If third-party items are used, functional failures may occur and fitness for use may be restricted. Biocompatibility requirements may also not be met. Please note that in such cases, any claim under warranty and liability will be voided if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- Use only nasal masks that do not have their own exhalation openings.
- Do not use masks that cover the entire nose and mouth.
- The muffling system has to be directly connected to the mask.
- Make sure that the connections between muffling system and the mask and between muffling system and the breathing hose are firm and air-tight. The connections must not be allowed to loosen during therapy.
- Do not close any openings on muffling system. Do not modify porous filter **2** in any way.
- If resistance should increase when exhaling, in spite of the muffling system having been cleaned regularly, replace the muffling system.
- Place the nasal mask on your face for any length of time only when the therapy device has been switched on. One of the purposes of the therapy device air flow is to remove the exhaled air from the mask and hose. When the therapy device is not switched on, there is a risk of rebreathing the exhaled air if you breathe through your nose.
- Do not use any electrically conductive breathing hoses.
- Clean the exhalation system as described in Point 4 "Hygienic preparation". Unsuitable cleaning methods, unsuitable detergents and insufficient cleaning frequency can lead to infection or bacterial contamination and impair the function of the exhalation system.

3 Operation

Putting into operation

Plug together the mask, the muffling system and the air hose of the treatment equipment.

After use

Disconnect the muffling system from the mask and the air hose. Clean the muffling system as described in the section entitled 4 "Hygienic preparation".

Combination with other equipment

The muffling system can be combined with CPAP and Bi-Level equipment when the following requirements are fulfilled:

- When the muffling system is combined with other devices, ensure that the therapy pressure and the pressure constancy in the nasal mask are as prescribed by your physician. Ask your physician or dealer whether the correct therapy pressure will be achieved in the mask using the desired equipment combination. Have your physician or dealer check the settings of your therapy device.
- Follow the directions in the user's manuals of the devices you want to combine with the muffling system.
- Use a nasal mask with a tapered connection in accordance with EN ISO 5356-1, e.g., SOMNOmask.
- Use a breathing hose made of elastomer with a connection diameter of 20 to 22 mm.

Cold/hot humidifiers

The muffling system can be combined with a heated or an unheated humidifier. Ensure that the muffling system is dry before starting treatment. Never select a degree of humidification which is so high that water forms in the breathing hose or exhalation system. Refer to the user's manuals for your humidifier and therapy device.

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4 Hygienic preparation

Intervals

A contaminated muffling system makes exhaling more difficult, so it is extremely important to clean the muffling system to ensure that the mask works properly.

Please clean the muffling system **after every use**.

Procedures

Disconnect the muffling system from the mask and the air hose. Clean the complete muffling system in clean warm water. Please do **not** use soap to clean the muffling system as it would stick the fine pores of the grey filter **2** together.

Dry the muffling system thoroughly after cleaning. Connect the air hose to your treatment equipment. Position the muffling system on the free end of the air hose and seal the open end of the muffling system with red plug **4**. Switch on the equipment and allow it to run for at least 30 minutes to ensure that it is thoroughly dry. After this time, switch off the equipment and **remove red plug 4**.

Change of patient

Use a new muffling system for every new patient.

5 Maintenance

The muffling system is maintenance-free.

Disposal

Dispose of used muffling systems via hospital or domestic waste disposal.

6 Scope of supply

Muffling system	WM 23685
Consisting of:	
– muffling system	WM 23625
– plug	WM 23639
– operating instructions	WM 16243

7 Technical data


Product category according to 93/42/EEC	IIa
Connections	
– Connecting piece	to fit connection with 22 mm Ø conical connector conforming to EN ISO 5356-1
– Sleeve	22.7 mm diameter, 1° external taper, suitable for elastomeric hoses with a connection diameter of 20 to 22 mm
Treatment pressure range	4 - 40 hPa
Weight	20 g
Temperature ranges:	
– Operation	+5 °C to +40 °C
– Storage	–20 °C to +70 °C
Max. moisture in hose	95 % rH not condensed
Flow resistance	
– at 50 l/min	0,07 hPa
– at 100 l/min	0,25 hPa
Average sound pressure level at 10 hPa	23,4 dB(A)
Duration of use	3 months, with careful maintenance up to 6 months, depending on the frequency of use and cleaning
Materials	
– Connecting piece	PC (polycarbonate)
– Sleeve	PC (polycarbonate)
– Porous filter	PE (polyethylene)

Subject to technical change without notice

Pressure/flow characteristic curve

The pressure/flow characteristic curve shows the outlet flow as a function of treatment pressure.

Muffling system average leakage flow [l/min.]

8 Warranty

Löwenstein Medical gives the customer a limited manufacturer warranty on new genuine Löwenstein Medical products and any replacement part fitted by Löwenstein Medical in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase as listed below. The warranty conditions can be downloaded from www.loewensteinmedical.de on the Internet. We can also send you the warranty conditions on request. In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty period
Masks including accessories, rechargeable battery, batteries (unless there is information to the contrary in the technical documentation), sensors, tube systems	6 months
Disposable products	None

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Technology GmbH & Co. KG