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Hersteller

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⚠ WARNING

The reconditioning of single use products can negatively influence the product properties. A change in the surface structure, for example, can lead to formation of cracks or a reduction in the tear strength as well as to microbiological aggregation of e.g. spores, allergens or pyrogens. Changing the material properties can lead to release of particles.

Risk of damage to patient's health!

- Only use single use products once and dispose of them properly after first use.

⚠ WARNING

Defective or leaking expiratory valves may not be used.

Risk of damage to patient's health!

- Visually inspect the expiratory valve and the diaphragm before use to ensure that it is undamaged.
- Dispose of the expiratory valve in case of damage or if the functional/leak test fails.

⚠ WARNING

Modifications to the product.

Risk of damage to patient's health!

- All modifications to the expiratory valve are prohibited.

NOTICE

- Follow the instructions for use of the ventilator.
- The responsible organization is liable for ensuring the compatibility of the ventilator with all other parts that are used for connecting to the patients before use.
- When using a plastic expiratory valve for the first time, the ventilator must be calibrated to these valves. The calibration may only be performed by a qualified technician. Please contact Service!
- Replace the expiratory valve as necessary according to your hospital's infection control procedures or in light of secretion deposits or atomization of medication.
- Löwenstein medical will provide no guarantee for the proper functioning of the single use expiratory valve if it is re-used.
- Every incident with the ventilator which leads to severe injury or death of the patient or to a potential public health hazard must be reported to the manufacturer and to the authorities responsible.

Intended use

The single use expiratory valve is intended for single use in Löwenstein medical ventilators. It serves to control the expiratory gas flow and maintain the pressure level set by the user.

To be noted for nebulization/humidification

If filters that are unsuitable for nebulization are used in the expiration side, they may become blocked, resulting in a significant increase in the flow resistance and ventilation hindrance. On the other hand, dispensing with the filters in the case of clouding with higher NaCl concentrations (e.g. 3%) results in increased salt deposits in the hose system and in the expiratory valve.

- Observe the instructions for use of the ventilator for using filters for nebulization.

Inspection prior to use

- Run through the function/leakage test according to the instructions for use of the ventilator after attaching the expiratory valve.

Fitting the single use expiratory valve

- Insert the diaphragm (1) in the diaphragm holder (2). The upper side of the diaphragm marked "TOP" must face outwards.

- Insert the expiratory valve in the correct position in the front of the appliance and turn it clockwise to the stop.

Scope of supply

The single use expiratory valve is delivered with the diaphragm clean and ready for clinical use.

Technical data

Conditions in operation: Temperature -10°C – 45°C; relative humidity 10% – 95%, no condensation.

Conditions in storage/transport: Temperature -20°C – 60°C; relative humidity 10% – 95%, no condensation.

Conical connectors: AD22/ID15 according to ISO 5356-1.

Compatible ventilators

- LeoniPlus / LeoniPlus transport
- Leoni 2 / Leoni m
- Leoni 4

Disposal

Used expiratory valves can be contaminated. Dispose of the expiratory valve according to your hospital's regulations and infection control procedures as well as the valid, local waste disposal regulations.

Symbols

Do not use if the packaging is damaged

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© 2023 Löwenstein Medical SE & Co. KG | Printed in Germany | All rights reserved
Löwenstein Medical Schweiz AG