

PesoCath

The key to esophageal and
transpulmonary pressure monitoring



Lung-protective through continuous Peso measurements

Measuring the TPP supplies up-to-date information about the mechanical pressure and volume load of the lungs and enables the ventilation regimen to be consistently adjusted for optimal lung protection.

Talmor D, Fessler H (2010) Are Esophageal Pressure Measurements Important in Clinical Decision-Making in Mechanically Ventilated Patients? Respiratory Care 55: 162-172

Monitoring esophageal pressure aids in identifying and treating causes of ineffectual patient efforts.

Beck J, Sinderby C, Lindström L, Grassino A (1998) Crural diaphragm activation during dynamic contractions at various inspiratory flow rates. J Appl Physiol 85: 451-458

Peso measurements of the work of breathing to quantify the respiratory effort enable the degree of muscle relief under ventilation to be individually adjusted to the patient's needs.

Brochard L (2014) Measurement of esophageal pressure at bedside: pros and cons. Curr Opin Crit Care 20: 39-46

In contrast to PEEPi measurements using a special ventilator manoeuvre, the intrinsic PEEP can be continuously measured in a spontaneously breathing patient with a correctly placed Peso catheter.


Zakynthinos SG, Vassilakopoulos T, Zakynthinos e, et al. (1997) Accurate measurements of intrinsic positive endexpiratory pressure: how to detect and correct for expiratory muscle activity. Eur Respir J 10:522-529





Wire-in-wire system

Easy withdrawal of guide wire through wire-in-wire solution and special coating



Esophageal balloon

The special esophageal balloon enables excellent response to erratic changes in pressure and is designed for the dynamic requirements of transpulmonary pressure measurement



Ventilator interface

Connection for esophageal and transpulmonary pressure monitoring or, depending on position, for gastric pressure



Simply safe in emergency

Direct channel for emergency suctioning, for auscultation and for draining gastric juices




Enteral nutrition, latest standard


Direct connection to the new standard for connectors according to DIN EN ISO 80369 with practical crown cap

What's what?

PesoCath esophageal catheter with:


- 1** Measurement channel
- 2** Guide wire (wire-in-wire technology)
- 3** Working channel for inserting the guide wire, for auscultatory confirmation of the position in the stomach and for draining gastric juices
- 4** Feeding channel with ENFit connector
- 5** Balloon with centred non-metallic X-ray marker






Accessories

Adaptors for additional pressure measurement ports of various intensive care ventilators



Pressure measurement line with three-way stopcock and luer-lock connectors



Syringe for balloon inflation



Technical specifications

Technical Data

Order numbers	1 piece: AZ-380230-1 5 pieces: AZ-380230-5
Exclusive distribution by	Löwenstein Medical Innovation
CE	0482
Component	Material
Gastric tube	PUR
Accessories	PVC, ABS, PE
Guide wire	Medical-grade stainless steel, PTFE-coated
Connectors	ABS, PVC, PE
Latex-free	Latex-free
Biocompatibility	ISO 10993-1
Classification	Class IIa
Application	Single patient use, short-term use for 30 days
Packaging	DIN EN ISO 11607-1, -2
Sterile packaging, primary packaging	Blister, PA/PE, medical grade paper




Dim	30 cm
Intended use	
Stomach tube for nu	
Technical data	
Order number	
Exclusive distribution by	
CE	
Component	M
Gastric tube	PUR
Accessories	PVC, A
Guide wire	Medical
Connectors	PTFE-coat
Latex-free	ABS, PVC, P
Biocompatibility	Latex-free
Classification	ISO 10993-1
Application	Class IIa
	Single patient use, short-term use for 30 d

Technical Data Gastric Tube Pesoc

Esophageal catheter PesoCath

1 piece: AZ-380230-1
5 pieces: AZ-380230-5

Illustration



- 1 Measurement channel
- 2 Working channel
- 3 Feeding channel
- 4 Wire-in-wire retrieval system
- 5 Balloon with radiopaque markers

Dimensions

25 cm

3 cm

150 cm

63 cm

29 cm

nutrition and oesophageal pressure measurement

1 piece: AZ-380230-1
5 pieces: AZ-380230-5

Löwenstein Medical Innovation
0482

Material

ABS, PE
grade stainless steel,
E

lays



made in
GERMANY

LÖWENSTEIN


medical

Sales + Service

Löwenstein Medical
Arzbacher Straße 80
56130 Bad Ems, Germany
T. +49 2603 9600-0
F. +49 2603 9600-50
info@loewensteinmedical.com
loewensteinmedical.com



Medicoplast International GmbH
Heusweilerstraße 100
66557 Illingen, Germany
T. +49 6825 9424-0

 With people in mind



p10390en2210

© Subject to copyright protection.
Reproduction of any kind requires the explicit permission of Löwenstein Medical.

CE 0482