

Confirmation

We, the manufacturer Weinmann Geräte für Medizin GmbH + Co. KG, Kronsaalsweg 40, 22525 Hamburg, confirm that our product

ACCUVAC Basic, Home Care Suction Pump

is suitable for use in aircraft – under the condition that the cabin air pressure is not below 0.7 bar (10.5 pounds per square inch). This has been confirmed to us by the manufacturer of the rechargeable battery pack.

The ACCUVAC Basic fulfills the standards involving electromagnetic compatibility, EN 55011 (Radio Interference) and EN 61000-4-2 to 5 and 11 (Immunity to Interference) as well as the EG standard 95/54/EG (e1-mark).

Hamburg, February 2006

WEINMANN
Geräte für Medizin GmbH & Co. KG



Dipl.-Kfm. Mark Brandmeier
Product Manager Home Mechanical Ventilation & Oxygen Medicine