

Leistungsübersicht

Samoa lite, Samoa, Scala, Sonata



Funktion	Samoa lite Screening	Samoa Polygraphie	Scala Polygraphie	Sonata Polysomnographie
EEG	-	-	✓ Mittels ExG Kanäle	✓
EOG	-	-	✓ Mittels ExG Kanäle	✓
EMG Kinn	-	-	✓ Mittels ExG Kanäle	✓
EKG	-	-	✓ Mittels ExG Kanäle	✓
LEMG	-	-	✓ Mittels ExG Kanäle	✓
ExG*	-	-	✓	✓
Flow (Nasenbrille)	✓	✓	✓	✓
Flow (Thermistor)	-	✓	✓	✓
Thorax/Abdomen (Druckpads)	-	✓	✓	✓
Thorax/Abdomen (RIP)	-	-	✓	✓
Schnarchen (Nasenbrille)	✓	✓	✓	✓
Schnarchen (Mikrofon)	-	-	✓	✓
SpO ₂	✓	✓	✓	✓
Puls	✓	✓	✓	✓
Pulswelle	✓	✓	✓	✓
PTT	-	-	-	✓
Lage	✓	✓	✓	✓
Licht	✓	✓	✓	✓
Druck	✓	✓	✓	✓
Atemfrequenz	✓	✓	✓	✓
Herzfrequenz	-	-	✓	✓
Phase T-A	-	✓	✓	✓
Blutdruck	-	-	-	✓
Aktimeter (externer Sensor)	-	-	-	✓
Aktimeter (intern)	✓	✓	-	-
Offline Messungen	✓	✓	✓	✓
Online Messungen	✓	✓	✓	✓
Digitale prisma Kanäle (online)	✓	✓	✓	✓
Analoge externe Kanäle	✓	✓	✓	✓
Audio / Video (online)	✓	✓	✓	✓
Bluetooth Kopplung (offline) zu prisma SMART / SOFT max / plus	-	✓	-	-
WLAN-fähig (online)	-	-	✓	✓

*2 Kanäle zur freien Konfiguration verfügbar z.B. EEG, EOG, EMG, LEMG oder EKG möglich.