

**LÖWENSTEIN**  
medical



# Allegra M30



The new generation of devices for secretion suction.

# Allegra M30

The Allegra M30 is future-oriented thanks to its modular design. The mechanical separation of the drive unit from the vacuum generator makes the Allegra M30 easy to service.

It combines technical innovation with user-oriented handling, giving users a new sense of security.

- Modern, attractive design
- High standard of hygiene: A patient-specific vacuum generator and an additional bacteria filter offer double safety
- Electronic monitoring of the suction force
- Simple operating concept

## Our highlights

- Innovative, modular design
- Container adapted directly to the device
- Mechanical separation of drive unit and vacuum generator
- Double safety thanks to additional bacteria filter
- Reliable locking of the container lid



## Unique reprocessing concept

In addition to its modern design and ease of use, the Allegra M30 offers a unique concept in which the modular construction mechanically separates the drive unit from the vacuum pump.

The vacuum generator is additionally protected by the overflow protection integrated in the container and the bacteria filter, thus enabling simple and cost-efficient reprocessing.

In addition to the reusable container, there is also a disposable container set, which enables user-friendly, quick suction and disposal.



### Technical specifications

Unit suction capacity (measuring point hose nozzle)	max. 30 l/min $\pm$ 2 l/
Vacuum	max. -90 kPa
Container	Reusable secretion canister system (1 l) Disposable secretion canister system (1 l)
Suction hose	CH25, length 150 cm (incl. fingertip)
CE marking	CE0494
Basic device dimensions (H x W x D)	213 x 260 x 124 mm
Weight (basic device)	2.8 kg
Operating time	5 minutes ON, 5 minutes OFF
Item number	REF 1510220h1 (incl. disposable container system) REF 1510215h1 (incl. reusable container system)



# LÖWENSTEIN medical

## Sales

Löwenstein Medical  
Arzbacher Straße 80,  
56130 Bad Ems, Germany  
T. +49 2603 9600-0  
F. +49 2603 9600-50  
info@loewensteinmedical.com  
loewensteinmedical.com



Asskea  
Haßlocher Str. 9  
99189 Gebesee, Germany  
T. +49 36201 57970  
info@asskea.de  
asskea.de



p10117en2504

© Protected by copyright. Reproduction of any kind requires the express permission of Löwenstein Medical.  
All information is subject to change for technical reasons.